



Instructions for Use
for
XOresearch Cardio.AI™

Software Version: 2.5

Date of Document Creation: 07-09-2023

Last Updated: 11-07-2025

Version: 1.5



Manufacturer:

XOresearch SIA

Republikas Square 3-225, Office 107, Riga, LV-1010, Latvia

Contact Information:

Phone: +371-67-305-084

Email: getintouch@xoresearch.com

Riga, Latvia
2025

Instruction for Use (IFU)

for XOresearch Cardio.AI™

Date	Version	Status/revision
07/09/2023	1.0	Document Creation
13/06/2024	1.1	Document update: Smart Actions section added to the document, instructions text has been adjusted.
23/09/2024	1.2	Document update: opening checklist removed, adjusted instruction text after internal reviewing.
28/03/2025	1.3	Document update: added sections "Availability of the Instructions for Use", "Limitations", legal address of the manufacturer clarified. Document update: added "Latvia", changed the year to 2025 under the first page. New limitation added. Software intended use has been changed IFU symbol has been changed to eIFU.
23/05/2025	1.4	Intended use change. Added sections "Awaiting For Record section overview". Added details of the expected lifetime.
11/07/2025	1.5	Performance characteristics have been formatted.

Manufacturer: XOresearch SIA.

Address: 3, Republikas square, office 107, Riga LV-1010, Latvia

Contact Information: Email: getintouch@xoresearch.com | Phone: +371-67-305-084








Software Identification:

- Software Name: XOresearch Cardio.AI™
- Software Version: Version 2.5
- Class IIa Software under Rule 11 of MDR (EU) 2017/745
- Basic UDI-DI: 47510473CARDIOAIQE | UDI-DI: 04751047370019



- **Intended Use:** XOresearch Cardio.AI™ is a standalone medical device software, using AI algorithms, intended to be used to analyze ECG signal recordings of adult patients by trained healthcare professionals in a clinical environment for the assessment of basic heart rate patterns. The results and interpretations generated by the software are reviewed, modified and approved by the physician. The physician retains full responsibility for the diagnosis and treatment decisions.

Symbols

	Manufacturer
 eIFU indicator	Consult electronic instructions for use (IFU)
	Medical device
	CE Mark and the Notified Body number
	Caution
	Catalogue number
	Unique Device Identifier



<div style="border: 2px solid black; padding: 5px; display: inline-block;"> SN </div>	Serial Number
--	---------------

Table of Contents:

Symbols	3
Table of Contents:	4
1. Introduction:	7
2. Software Description:	7
3. Indications for Use:	8
4. Contraindications/Warnings:	8
5. Patient Population	9
6. Intended users:	9
7. Precautions/Warnings:	9
8. Expected Lifetime	10
9. Performance characteristics of the device	10
9.1 Accuracy	11
9.2 AUC	14
9.3 F1-Score	14
9.4 PPV	17
9.5 Sensitivity	20
9.6 Specificity	23
9.7 False Negatives	26
9.8 False Positives	29
10. Technical Requirements:	31
11. Setup:	32
12. Software Operation:	33
12.1 User profile management	33
12.1.1 Edit user data	34
12.1.2 Update user data	36
12.1.3 Delegate control of my account to a user	37
12.1.4 User interface language	39
12.1.5 Hiding sensitive information	40
12.2 Organizations overview	41
12.2.1 Organizations filters	42
12.2.2 Calculate spendings of the organization	43



12.2.3 Export spendings into CSV	43
12.2.4 Edit organization	44
12.2.5 Organization report preset configuration	47
12.2.6 Remove organization	52
12.3 Awaiting For Record section overview	53
12.3.1 Messenger access	55
12.3.2 Binding device to patient	58
12.3.3 Patient data editing	61
12.4 Tasks section overview	61
12.4.1 Reviewing subsection overview	61
12.4.2 Reviewing subsection editing	64
12.4.2.1 Patient data editing	64
12.4.2.2 Channels editing	66
12.4.2.3 Task reclassification	72
12.4.2.3 Replacing task data	72
12.4.2.4 Task cancellation	73
12.4.2.5 Task deletion	73
12.4.3 Uploading subsection overview	74
12.5 Users section	76
12.5.1 Users section overview	76
12.5.2 User creation	78
12.5.3 User invitation	80
12.5.4 User editing	81
12.5.5 User role assignment deletion	81
12.6 Roles section	82
12.6.1 Roles section overview	82
12.6.2 Role management	86
12.7 ECG data input	88
12.8 ECG data analysing	89
12.8.1 ECG Viewer Header	90
12.8.1.1 Share ECG task	91
12.8.1.2 ECG task options	91
12.8.1.3 Edit ECG task	95
12.8.1.4 Close ECG task	95
12.8.1.5 Pre-approve ECG task	96
12.8.1.6 Approve ECG task	96
12.8.1.7 Save ECG task	96
12.8.2 ECG Viewer Editor	96
12.8.2.1 Side editing bar	97
12.8.2.2 Poincare plot	99

12.8.2.3 Beats list	103
12.8.2.4 Smart actions	106
12.8.2.5 Beats clusters panel	109
12.8.2.6 Beats cross annotations list	110
12.5.3 ECG Viewer Previewer	111
12.8.4 ECG Viewer Visualizer	112
12.5.5 ECG Viewer bird view	124
12.9 ECG data report	124
12.9.1 ECG data report overview	124
12.9.2 ECG data report sections managing	132
12.9.3 ECG data report editing	134
13. Data Input and Output:	135
14. User Authentication and Access Control:	135
15. Data Security and Privacy:	137
16. Troubleshooting:	138
17. Availability of the Instructions for Use (IFU):	138
18. Limitations	139
19. Manufacturer's Declaration	141



1. Introduction:

Welcome to the Instruction for Use (IFU) for XOresearch Cardio.AI™. This document is provided by XOresearch SIA to assist healthcare professionals in safely and effectively utilizing our clinical decision support software.

The IFU contains essential information about the software's features, its intended use, precautions, and guidance on troubleshooting. Please read this document carefully before using the software.

2. Software Description:

XOresearch Cardio.AI™ is a multipurpose device for automatic annotation and interpretation primarily long, and short (from 7 seconds to 35 days) ECG records with any lead combinations and designed to:

- detect heartbeats, in pre-recorded ECG data; and separate noise from the beats in the data analysed by the device; and
- detect beat and rhythm events for the following rhythms: sinus, atrial, junctional, ventricular; and for the following signs of: pre-excitation, conduction abnormalities, intraventricular conduction delays;; and
- detect PQRST points, ST segment amplitude and direction, T-wave type, HRV, Heart Rate BPM; and
- visualise ECG data along with the other vital signs and patient-related information such as indications, diary events, demographic data; and
- generate an interpretation statement on an ECG data; and
- create a report based on the ECG findings in export it in a PDF format alongside priority indicators labelling; and
- store ECG data in cloud storage; and
- provide temporary or permanent access to ECG data or other vital signs.

The annotation made by the device will be confirmed by the physician and can be edited, or deleted. Platform interpretation results are not intended to be the sole means of diagnosis.

XOresearch Cardio.AI™ is a multipurpose medical device designed by the manufacturer to serve the following clinical purposes:

Automatic Annotation and Interpretation: The primary function of this device is the automatic annotation and interpretation of primarily long and short ECG records, irrespective of lead combinations.

It is specifically developed to:

Detect Heartbeats: Accurately identify heartbeats within pre-recorded ECG data.



Noise Separation: Distinguish and separate noise from the analyzed heartbeats in the data.

Rhythm Detection: Detect various cardiac rhythms, including sinus, atrial, junctional, and ventricular rhythms.

Disorder Identification: Identify specific cardiac disorders such as pre-excitation syndromes, heart blocks, and bundle branch blocks.

Data Analysis: Analyze critical ECG parameters like PQRST points, ST segment amplitude and direction, T-wave type, Heart Rate Variability (HRV), and Heart Rate in beats per minute (BPM).

Comprehensive Visualization: Display ECG data alongside vital signs and patient-related information, including indications, diary events, and demographic data.

Interpretation Generation: Generate an interpretation statement based on the analyzed ECG data.

Report Creation: Create a comprehensive report summarizing ECG findings, exportable in PDF format, with severity indicators labeling.

Cloud Storage: Store ECG data securely in cloud storage for easy access and retrieval.

Data Accessibility: Provide both temporary and permanent access to ECG data and other vital signs as required.

Please note that while the device offers automatic annotation and interpretation, it is crucial to emphasize that these results are not intended to serve as the sole means of diagnosis. Physicians may confirm, edit, or delete annotations made by the device as part of their clinical practice.

XOresearch Cardio.AI™ software supports import of EDF and BDF data files from compatible ECG Holter devices via manual upload and API-based transfer methods.

3. Indications for Use:

- XOresearch Cardio.AI™ is intended for use in a hospital or clinical setting, by or on the order of a physician or similarly qualified health care professional. XOresearch Cardio.AI™ evaluates the ECG data of ambulatory patients pre-recorded with a legally marketed digital ECG recorder with any lead combinations.

The annotation by the device will be confirmed and may be edited or deleted by the physician. The final decision regarding the treatment of patients is the responsibility of the physician. Platform interpretation results are not intended to be the sole means of diagnosis.

4. Contraindications/Warnings:

XOresearch Cardio.AI™ is not indicated to detect pacemaker, as pacemaker's detection is not part of the current version of the system. XOresearch Cardio.AI™ does not analyse the pacemaker function and threats the signal as is, without any assumptions on presence or absence of the pacemaker, thus should not be used in fully automatic mode without physician attention for patients with a pacemaker.



XOresearch Cardio.AI™ does not support online real time analysis of ECG data. XOresearch Cardio.AI™ processes offline data in post-processing mode.
XOresearch Cardio.AI™ is not intended for real-time patient monitoring.

5. Patient Population

XOresearch Cardio.AI™ intended to be used on adult patient records (over the age of eighteen) prescribed to undergo electrocardiography.

6. Intended users:

XOresearch Cardio.AI™ is intended for use by medical professionals, such as those who are responsible for deciphering ECG data, analyzing it, and diagnosing the patient on the basis of this data.

Operators of XOresearch Cardio.AI™ must hold recognized qualifications in cardiology or a related discipline, per Directive 2005/36/EC.

All operators of XOresearch Cardio.AI™ should thoroughly read and acknowledge this IFU to ensure safe and effective use. Acknowledgment of the IFU confirms that the user understands the capabilities, limitations, and best practices associated with the software.

7. Precautions/Warnings:

- Ensure that your computer system meets the minimum system requirements specified in the software's documentation. Inadequate hardware or software configurations may result in performance issues or software malfunctions.
- Verify the accuracy of input data, as inaccurate or incomplete data may lead to incorrect recommendations.
- Use the software in a controlled clinical environment with proper lighting and minimal distractions to minimize the risk of errors.
- Always exercise sound clinical judgment when interpreting the software's recommendations. The software is a decision support tool and should not replace the expertise of healthcare professionals.
- Do not rely solely on the software's recommendations for critical or life-threatening decisions. In such cases, seek immediate clinical assessment and intervention.
- Report any software-related issues, errors, or discrepancies to the appropriate personnel or IT support to address and resolve them promptly.
- Ensure that healthcare professionals using the software are adequately trained and competent in its use. Training should cover software operation, data input, result interpretation, and troubleshooting.
- Do not solely rely on the software's recommendations; use clinical judgment.
- Ensure data input is accurate, as incorrect data may lead to incorrect recommendations.



- XOresearch Cardio.AITM is a decision support tool and is not a substitute for the expertise of trained healthcare professionals. Healthcare providers must exercise their clinical judgment when interpreting software recommendations and making medical decisions.
- In cases of urgent or critical medical conditions where immediate clinical assessment and intervention are required, do not solely rely on the software's recommendations. Delaying necessary actions can have serious consequences.
- The accuracy of outputs generated by the software depends on the accuracy and completeness of the input data. Users are responsible for verifying the correctness of patient data entered into the system.
- Healthcare professionals are solely responsible for interpreting and acting upon the software's recommendations. Exercise caution and ensure that the recommendations align with the patient's clinical presentation and history.
- Protect patient data and ensure its security during transmission and storage. Unauthorized access or data breaches can compromise patient privacy and confidentiality.
- Report any software-related errors, discrepancies, or unusual behaviors to your organization's IT support or the software provider. Do not attempt to modify or alter the software without proper authorization.
- The software provider and manufacturer disclaim liability for any adverse events or consequences arising from the use of XOresearch Cardio.AITM beyond the extent permitted by law. Healthcare professionals are responsible for their decisions and actions.

8. Expected Lifetime

XOresearch Cardio.AI™ is designed to maintain its intended performance and safety for an expected lifetime of 15 years, under normal operating conditions and with proper maintenance. The lifetime includes continued support through software updates, cybersecurity patches, and compatibility adaptations, as necessary to meet evolving regulatory and technical standards.

9. Performance characteristics of the device

XOresearch Cardio.AI™ has the following performance metrics:

- Accuracy
- Area Under the Curve (AUC)
- F1-Score
- Positive Predictive Value (PPV)
- Sensitivity
- False Negatives
- False Positives

The device is expected to remain safe and effective for a period of 15 years, as specified by the manufacturer. This is based on lifecycle validation activities, including risk management,



state-of-the-art assessment, and post-market surveillance planning. Throughout this lifetime, XOresearch SIA commits to maintain the device's clinical performance through validated software updates.

9.1 Accuracy

Accuracy indicates the overall performance of the classification model by calculating the proportion of correctly predicted instances (both positives and negatives) out of the total number of instances. I

Label	Accuracy
Atrial Premature Contraction	0,9999
Aberrated Atrial Premature Beat	0,9999
Non-Conducted P-Wave (Blocked)	0,9975
Left Anterior Fascicular Block Beat (Common)	0,9999
Bifascicular Block Beat	0,9999
Intraventricular Conduction Disturbance (Non-Specific Block)	0,9999
Left Posterior Fascicular Block Beat (Rare)	0,9999
Junctional (Nodal) Escape Beat	0,9993
Junctional (Nodal) Premature Contraction	0,9782
Left Bundle Branch Block Beat	0,9999
Incomplete Left Bundle Branch Block Beat	0,9999
Normal Beat	0,9999
Right Bundle Branch Block Beat	0,9999
Incomplete Right Bundle Branch Block Beat	0,9999



Unclassifiable Beat	0,9992
Ventricular Escape Beat	0,9896
Ventricular Premature Contraction	0,9999
Fusion Of Ventricular And Normal Beat	0,9979
Noise (No Signal)	0,9999
Noise Severe	0,9999
Asystole	1
Atrial Ectopic Rhythm	0,9999
Atrial Fibrillation	0,9999
Atrial Flutter	0,9999
Multifocal Atrial Tachycardia	0,9999
Paroxysmal Atrial Tachycardia	0,9999
AV Dissociation With Interference	1
First Degree AV Block	0,9999
Second Degree AV Block Type I	0,9993
Second Degree AV Block Type II	0,9999
Third Degree AV Block	0,9999
Accelerated Av Junctional (Nodal) Rhythm	0,9999
AV Junctional (Nodal) Escape Rhythm	0,9999
Junctional Tachycardia	0,9999
Lown-Ganong-Levine Syndrome	0,9999
Second Degree SA Block Type I	0,9998
Second Degree SA Block Type II	0,9995
Third Degree SA Block	0,9969



Sinus Arrhythmia	0,9999
Sinus Tachycardia	0,9999
Accelerated Idioventricular Rhythm	0,9996
Ventricular Fibrillation	0,9729
Idioventricular (Ventricular Escape) Rhythm	0,9913
Ventricular Couplet	0,9999
Monomorphic Ventricular Tachycardia	0,9999
Polymorphic Ventricular Tachycardia	0,9999
Torsades De Pointes Ventricular Tachycardia	0,9987
Wandering Pacemaker From The Sinus Node To (And From) The A-V Node	0,9999
Wolf-Parkinson Type A	1
Wolf-Parkinson Type B	0,9999
Artifact	0,9996
Ventricular Interpolated Beat	0,9991
Atrial Couplet	0,9999
Atrial Triplet	0,9999
Junctional Couplet	0,9945
Junctional Triplet	0,9999
Ventricular Triplet	0,9999



9.2 AUC

AUC (Area Under the Curve) represents the area under the Receiver Operating Characteristic (ROC) curve, which plots the True Positive Rate (Sensitivity) against the False Positive Rate (1-Specificity) at various threshold levels. AUC measures the model's ability to distinguish between positive and negative classes.

AUC value is **0.9991412278967556**

9.3 F1-Score

F1 score a balanced measure of a classification model's performance. It is especially useful when there is an uneven class distribution or when false positives and false negatives have different consequences.

Label	F1
Atrial Premature Contraction	0.9834
Aberrated Atrial Premature Beat	0.9634
Non-Conducted P-Wave (Blocked)	0.9512
Left Anterior Fascicular Block Beat (Common)	0.9999
Bifascicular Block Beat	0.8854
Intraventricular Conduction Disturbance (Non-Specific Block)	0.9986
Left Posterior Fascicular Block Beat (Rare)	0.9995
Junctional (Nodal) Escape Beat	0.939
Junctional (Nodal) Premature Contraction	0.7755
Left Bundle Branch Block Beat	0.9808
Incomplete Left Bundle Branch Block Beat	0.9992
Normal Beat	0.9975



Right Bundle Branch Block Beat	0.8914
Incomplete Right Bundle Branch Block Beat	0.9655
Unclassifiable Beat	0.9419
Ventricular Escape Beat	0.9143
Ventricular Premature Contraction	0.9923
Fusion Of Ventricular And Normal Beat	0.9189
Noise (No Signal)	0.9941
Noise Severe	0.9348
Asystole	1
Atrial Ectopic Rhythm	0.9948
Atrial Fibrillation	0.9996
Atrial Flutter	0.9818
Multifocal Atrial Tachycardia	0.959
Paroxysmal Atrial Tachycardia	0.9504
AV Dissociation With Interference	1
First Degree AV Block	0.9941
Second Degree AV Block Type I	0.9244
Second Degree AV Block Type II	0.9846
Third Degree AV Block	0.9965
Accelerated Av Junctional (Nodal) Rhythm	0.9964
AV Junctional (Nodal) Escape Rhythm	0.9924
Junctional Tachycardia	0.9799
Lown-Ganong-Levine Syndrome	0.9878
Second Degree SA Block Type I	0.9787



Second Degree SA Block Type II	0.968
Third Degree SA Block	0.9
Sinus Arrhythmia	0.9502
Sinus Tachycardia	0.9905
Accelerated Idioventricular Rhythm	0.9716
Ventricular Fibrillation	0.8571
Idioventricular (Ventricular Escape) Rhythm	0.9231
Ventricular Couplet	0.9936
Monomorphic Ventricular Tachycardia	0.9958
Polymorphic Ventricular Tachycardia	0.9248
Torsades De Pointes Ventricular Tachycardia	0.7481
Wandering Pacemaker From The Sinus Node To (And From) The A-V Node	0.9882
Wolf-Parkinson Type A	1
Wolf-Parkinson Type B	0.9985
Artifact	0.9623
Ventricular Interpolated Beat	0.9792
Atrial Couplet	0.9907
Atrial Triplet	0.9871
Junctional Couplet	0.8889
Junctional Triplet	0.9913
Ventricular Triplet	0.9857



9.4 PPV

Positive Predictive Value (PPV) represents the proportion of true positive predictions out of all instances that the model classified as positive.

Label	Precision
Atrial Premature Contraction	0.9754
Aberrated Atrial Premature Beat	0.9527
Non-Conducted P-Wave (Blocked)	1
Left Anterior Fascicular Block Beat (Common)	0.9999
Bifascicular Block Beat	0.7946
Intraventricular Conduction Disturbance (Non-Specific Block)	0.9982
Left Posterior Fascicular Block Beat (Rare)	0.999
Junctional (Nodal) Escape Beat	0.9365
Junctional (Nodal) Premature Contraction	0.9048
Left Bundle Branch Block Beat	0.9625
Incomplete Left Bundle Branch Block Beat	0.9996
Normal Beat	0.9981
Right Bundle Branch Block Beat	0.8045
Incomplete Right Bundle Branch Block Beat	1
Unclassifiable Beat	0.9625
Ventricular Escape Beat	0.9412
Ventricular Premature Contraction	0.9977



Fusion Of Ventricular And Normal Beat	0.8947
Noise (No Signal)	0.9912
Noise Severe	0.9275
Asystole	1
Atrial Ectopic Rhythm	0.9929
Atrial Fibrillation	0.9996
Atrial Flutter	0.9646
Multifocal Atrial Tachycardia	0.9915
Paroxysmal Atrial Tachycardia	0.9989
AV Dissociation With Interference	1
First Degree AV Block	0.9901
Second Degree AV Block Type I	0.9554
Second Degree AV Block Type II	0.9811
Third Degree AV Block	1
Accelerated Av Junctional (Nodal) Rhythm	0.9976
AV Junctional (Nodal) Escape Rhythm	1,0
Junctional Tachycardia	0.9841
Lown-Ganong-Levine Syndrome	0.9793
Second Degree SA Block Type I	0.9871
Second Degree SA Block Type II	1
Third Degree SA Block	0.9
Sinus Arrhythmia	0.9627
Sinus Tachycardia	0.9836
Accelerated Idioventricular Rhythm	1
Ventricular Fibrillation	0.75



Idioventricular (Ventricular Escape) Rhythm	1
Ventricular Couplet	0.9882
Monomorphic Ventricular Tachycardia	0.9949
Polymorphic Ventricular Tachycardia	0.9295
Torsades De Pointes Ventricular Tachycardia	0.6898
Wandering Pacemaker From The Sinus Node To (And From) The A-V Node	0.993
Wolf-Parkinson Type A	1
Wolf-Parkinson Type B	0.9975
Artifact	0.9746
Ventricular Interpolated Beat	0.9792
Atrial Couplet	0.9938
Atrial Triplet	0.9894
Junctional Couplet	0.9091
Junctional Triplet	0.9956
Ventricular Triplet	0.9942



9.5 Sensitivity

Sensitivity measures the proportion of correctly identified positive instances out of all actual positive instances.

Label	Sensitivity
Atrial Premature Contraction	0.9916
Aberrated Atrial Premature Beat	0.9743
Non-Conducted P-Wave (Blocked)	0.907
Left Anterior Fascicular Block Beat (Common)	0.9999
Bifascicular Block Beat	0.9995
Intraventricular Conduction Disturbance (Non-Specific Block)	0.999
Left Posterior Fascicular Block Beat (Rare)	0.9999
Junctional (Nodal) Escape Beat	0.9415
Junctional (Nodal) Premature Contraction	0.6786
Left Bundle Branch Block Beat	0.9998
Incomplete Left Bundle Branch Block Beat	0.9988
Normal Beat	0.9969
Right Bundle Branch Block Beat	0.9993
Incomplete Right Bundle Branch Block Beat	0.9334
Unclassifiable Beat	0.9222
Ventricular Escape Beat	0.8889
Ventricular Premature Contraction	0.9869
Fusion Of Ventricular And Normal Beat	0.9444



Noise (No Signal)	0.9969
Noise Severe	0.9422
Asystole	1
Atrial Ectopic Rhythm	0.9967
Atrial Fibrillation	0.9997
Atrial Flutter	0.9996
Multifocal Atrial Tachycardia	0.9287
Paroxysmal Atrial Tachycardia	0.9064
AV Dissociation With Interference	1
First Degree AV Block	0.9982
Second Degree AV Block Type I	0.8954
Second Degree AV Block Type II	0.9882
Third Degree AV Block	0.993
Accelerated Av Junctional (Nodal) Rhythm	0.9952
AV Junctional (Nodal) Escape Rhythm	0.9849
Junctional Tachycardia	0.9757
Lown-Ganong-Levine Syndrome	0.9965
Second Degree SA Block Type I	0.9705
Second Degree SA Block Type II	0.9379
Third Degree SA Block	0.9
Sinus Arrhythmia	0.938
Sinus Tachycardia	0.9974
Accelerated Idioventricular Rhythm	0.9448
Ventricular Fibrillation	1
Idioventricular (Ventricular Escape)	0.8571



Rhythm	
Ventricular Couplet	0.999
Monomorphic Ventricular Tachycardia	0.9967
Polymorphic Ventricular Tachycardia	0.9201
Torsades De Pointes Ventricular Tachycardia	0.8172
Wandering Pacemaker From The Sinus Node To (And From) The A-V Node	0.9834
Wolf-Parkinson Type A	1
Wolf-Parkinson Type B	0.9995
Artifact	0.9504
Ventricular Interpolated Beat	0.9792
Atrial Couplet	0.9876
Atrial Triplet	0.9848
Junctional Couplet	0.8696
Junctional Triplet	0.987
Ventricular Triplet	0.9773



9.6 Specificity

Specificity measures the proportion of correctly identified negative instances out of all actual negative instances.

Label	Specificity
Atrial Premature Contraction	0,9991
Aberrated Atrial Premature Beat	0,9983
Non-Conducted P-Wave (Blocked)	1
Left Anterior Fascicular Block Beat (Common)	0
Bifascicular Block Beat	0,9643
Intraventricular Conduction Disturbance (Non-Specific Block)	0,9999
Left Posterior Fascicular Block Beat (Rare)	0,9989
Junctional (Nodal) Escape Beat	0,9954
Junctional (Nodal) Premature Contraction	0,9851
Left Bundle Branch Block Beat	0,9871
Incomplete Left Bundle Branch Block Beat	0
Normal Beat	0,9999
Right Bundle Branch Block Beat	0,9737
Incomplete Right Bundle Branch Block Beat	1
Unclassifiable Beat	0,9974
Ventricular Escape Beat	0,9697
Ventricular Premature Contraction	0,9999



Fusion Of Ventricular And Normal Beat	0,9622
Noise (No Signal)	0,9985
Noise Severe	0,9995
Asystole	N/A
Atrial Ectopic Rhythm	0,9989
Atrial Fibrillation	0,9999
Atrial Flutter	0,9909
Multifocal Atrial Tachycardia	0,9999
Paroxysmal Atrial Tachycardia	0,9999
AV Dissociation With Interference	N/A
First Degree AV Block	0,9997
Second Degree AV Block Type I	0,9976
Second Degree AV Block Type II	0,998
Third Degree AV Block	1
Accelerated Av Junctional (Nodal) Rhythm	0,9995
AV Junctional (Nodal) Escape Rhythm	1
Junctional Tachycardia	0,9986
Lown-Ganong-Levine Syndrome	0,9989
Second Degree SA Block Type I	0,9935
Second Degree SA Block Type II	1
Third Degree SA Block	0,9
Sinus Arrhythmia	0,9993
Sinus Tachycardia	0,9997
Accelerated Idioventricular Rhythm	1
Ventricular Fibrillation	0



Idioventricular (Ventricular Escape) Rhythm	1
Ventricular Couplet	0,9881
Monomorphic Ventricular Tachycardia	0,9994
Polymorphic Ventricular Tachycardia	0,9916
Torsades De Pointes Ventricular Tachycardia	0,9955
Wandering Pacemaker From The Sinus Node To (And From) The A-V Node	0,9999
Wolf-Parkinson Type A	N/A
Wolf-Parkinson Type B	0,9974
Artifact	0,9983
Ventricular Interpolated Beat	0,9792
Atrial Couplet	0,9997
Atrial Triplet	0,9988
Junctional Couplet	0,9836
Junctional Triplet	0,9985
Ventricular Triplet	0,9995



9.7 False Negatives

False Negatives (FN) identifies instances where the classification model incorrectly predicts a positive instance as negative.

Label	False Negatives
Atrial Premature Contraction	0,0084
Aberrated Atrial Premature Beat	0,0263
Non-Conducted P-Wave (Blocked)	0,1025
Left Anterior Fascicular Block Beat (Common)	0,0001
Bifascicular Block Beat	0,0005
Intraventricular Conduction Disturbance (Non-Specific Block)	0,001
Left Posterior Fascicular Block Beat (Rare)	0,0001
Junctional (Nodal) Escape Beat	0,0621
Junctional (Nodal) Premature Contraction	0,4735
Left Bundle Branch Block Beat	0,0002
Incomplete Left Bundle Branch Block Beat	0,0012
Normal Beat	0,0031
Right Bundle Branch Block Beat	0,0007
Incomplete Right Bundle Branch Block Beat	0,0713
Unclassifiable Beat	0,0843
Ventricular Escape Beat	0,1249
Ventricular Premature Contraction	0,0132



Fusion Of Ventricular And Normal Beat	0,0588
Noise (No Signal)	0,0031
Noise Severe	0,0613
Asystole	0
Atrial Ectopic Rhythm	0,0033
Atrial Fibrillation	0,0003
Atrial Flutter	0,0004
Multifocal Atrial Tachycardia	0,0767
Paroxysmal Atrial Tachycardia	0,1032
AV Dissociation With Interference	0
First Degree AV Block	0,0018
Second Degree AV Block Type I	0,1168
Second Degree AV Block Type II	0,0119
Third Degree AV Block	0,007
Accelerated Av Junctional (Nodal) Rhythm	0,0048
AV Junctional (Nodal) Escape Rhythm	0,0153
Junctional Tachycardia	0,0249
Lown-Ganong-Levine Syndrome	0,0035
Second Degree SA Block Type I	0,0303
Second Degree SA Block Type II	0,0662
Third Degree SA Block	0,1111
Sinus Arrhythmia	0,066
Sinus Tachycardia	0,0026
Accelerated Idioventricular Rhythm	0,0584
Ventricular Fibrillation	0



Idioventricular (Ventricular Escape) Rhythm	0,1667
Ventricular Couplet	0,001
Monomorphic Ventricular Tachycardia	0,0033
Polymorphic Ventricular Tachycardia	0,0868
Torsades De Pointes Ventricular Tachycardia	0,2236
Wandering Pacemaker From The Sinus Node To (And From) The A-V Node	0,0168
Wolf-Parkinson Type A	0
Wolf-Parkinson Type B	0,0005
Artifact	0,0521
Ventricular Interpolated Beat	0,0212
Atrial Couplet	0,0125
Atrial Triplet	0,0154
Junctional Couplet	0,1499
Junctional Triplet	0,0131
Ventricular Triplet	0,0232



9.8 False Positives

False Positives (FP) identifies instances where the classification model incorrectly predicts a negative instance as positive.

Label	False Positives
Atrial Premature Contraction	0,0252
Aberrated Atrial Premature Beat	0,0496
Non-Conducted P-Wave (Blocked)	0
Left Anterior Fascicular Block Beat (Common)	0,0001
Bifascicular Block Beat	0,2585
Intraventricular Conduction Disturbance (Non-Specific Block)	0,0018
Left Posterior Fascicular Block Beat (Rare)	0,001
Junctional (Nodal) Escape Beat	0,0678
Junctional (Nodal) Premature Contraction	0,1052
Left Bundle Branch Block Beat	0,0389
Incomplete Left Bundle Branch Block Beat	0,0004
Normal Beat	0,0019
Right Bundle Branch Block Beat	0,243
Incomplete Right Bundle Branch Block Beat	0
Unclassifiable Beat	0,0389
Ventricular Escape Beat	0,0624
Ventricular Premature Contraction	0,0023



Fusion Of Ventricular And Normal Beat	0,1176
Noise (No Signal)	0,0088
Noise Severe	0,0781
Asystole	0
Atrial Ectopic Rhythm	0,0071
Atrial Fibrillation	0,0004
Atrial Flutter	0,0367
Multifocal Atrial Tachycardia	0,0085
Paroxysmal Atrial Tachycardia	0,0011
AV Dissociation With Interference	0
First Degree AV Block	0,0099
Second Degree AV Block Type I	0,0466
Second Degree AV Block Type II	0,0192
Third Degree AV Block	0
Accelerated Av Junctional (Nodal) Rhythm	0,0024
AV Junctional (Nodal) Escape Rhythm	0
Junctional Tachycardia	0,0161
Lown-Ganong-Levine Syndrome	0,0211
Second Degree SA Block Type I	0,013
Second Degree SA Block Type II	0
Third Degree SA Block	0,1111
Sinus Arrhythmia	0,0387
Sinus Tachycardia	0,0166
Accelerated Idioventricular Rhythm	0
Ventricular Fibrillation	0,3333



Idioventricular (Ventricular Escape) Rhythm	0
Ventricular Couplet	0,0119
Monomorphic Ventricular Tachycardia	0,0051
Polymorphic Ventricular Tachycardia	0,0758
Torsades De Pointes Ventricular Tachycardia	0,4496
Wandering Pacemaker From The Sinus Node To (And From) The A-V Node	0,007
Wolf-Parkinson Type A	0
Wolf-Parkinson Type B	0,0025
Artifact	0,026
Ventricular Interpolated Beat	0,0212
Atrial Couplet	0,0062
Atrial Triplet	0,0107
Junctional Couplet	0,0999
Junctional Triplet	0,0044
Ventricular Triplet	0,0058

10. Technical Requirements:

XOresearch Cardio.AI™ is accessible through a web browser based on Chromium browser engine: Google Chrome, Microsoft Edge, Opera Browser.

It is recommended to use the latest version of Google Chrome for optimal compatibility and performance.

The minimum required of the last stable version Google Chrome is 116 - when the IFU was produced).

The minimum required of the last stable version of Microsoft Edge is 126, of the Opera - 113.

- A stable and high-speed internet connection is essential for accessing XOresearch Cardio.AI™. A minimum download and upload speed of 100Mbit/s is recommended.



- Ensure that network firewall and security settings allow access to the XOresearch Cardio.AI™ web application. It may be necessary to whitelist the following domains to ensure unobstructed access: <https://web.cardio.ai/>

Operating System: XOresearch Cardio.AI™ is compatible with Windows 11, 22H2, macOS .

Hardware minimal requirements for running Google Chrome to access XOresearch Cardio.AI™ are:

Processor: 1.6 GHz or faster processor (Intel Pentium 4 or later).

RAM: 2 GB (minimum) for normal usage, 4 GB or more recommended for better performance.

Hard Drive: At least 100 MB of free space for browser installation.


Graphics: Graphics hardware acceleration requires a DirectX 9.0c capable video card with WDDM 1.0 or higher driver.

The electronic Instructions for Use (eIFU) are hosted on a secure and reliable web-based platform with high availability. The server infrastructure ensures minimal downtime, and access to IFU remains uninterrupted under normal operating conditions. Users experiencing accessibility issues should contact technical support at getintouch@xoresearch.com.

11. Setup:


- The option to access XOresearch Cardio.AI™ is available under the following web link: <https://web.cardio.ai/>

Software shows the following screen when successful:




The sign in to XOresearch Cardio.AI™ is available under the filling in the **Email or phone number** field > **Password** field > Continue button:





Note: Login credentials are provided by the manufacturer.

XOresearch Cardio.AI™ shows the following screen when successful sign in:




Note. The user is automatically logged out after 10 minutes of inactivity.

12. Software Operation:


12.1 User profile management

The option to access User profile management is available by clicking on the user name:



12.1.1 Edit user data

The option to edit user data is available under User profile management menu > My settings:



XOresearch Cardio.AI™ shows the following screen when successful:



Update user data

First and Last name*

Demis Constantinos

Current Password*

New Password*



Confirm New Password*

Ordering Organization

Ordering Phone

Ordering Address

Cancel

Update

Delegate control of my account to a user

Email

Add manager

User token

Token

Token does not exist.

The following sections are presented under the **My settings**:

- Update user data;
- Delegate control of my account to a user;
- User token.



12.1.2 Update user data

The following settings are available to update under **Update user data** section (All required fields are marked with an asterisk *):

Setting	Description
First and Last name*	Indicates the First and Last name of the user, visible. This field is required .
Current password*	Enables to provide the current password in order to change it. This field is required when password changing.
New password*	Indicates user's password that will be used during the log in process. Password requirements: <ul style="list-style-type: none">• At least 1 special symbol;• At least 1 lowercase letter;• At least 1 uppercase letter;• At least 1 digit;• Length must be at least 8 symbols. This field is required when password changing.
Confirm new password*	This field duplicates the Password field and must be filled in identically. This field is required when password changing.
Ordering Organization	This field indicates the organization which the user is associated with.
Ordering Phone	This field indicates the phone number which the user is associated with.
Ordering Address	This field indicates the address which the user is associated with.
Notification email	This field indicates the email address the generated task report will be sent to.

The option to update the data is available by filling in the data to the correspondent field, and clicking **Update** button. The option to Cancel changes and close the windows is available under the **Cancel** button.

The option to update password is available by filling in the **Current password**, **New password** and **Confirm new password** fields, and clicking the **Update** button.




The image shows two versions of a password update form. The left version is the initial state with empty fields: 'Current Password*', 'New Password*' (with a radio button), 'Confirm New Password*', 'Ordering Organization', 'Ordering Phone', and 'Ordering Address'. Below are 'Cancel' and 'Update' buttons. The right version shows the form after a password change, with the password fields filled with dots. A green arrow points from the left form to the right form, indicating a transition or update.

12.1.3 Delegate control of my account to a user

Delegating control of the account to another user enables a user to operate the delegated account by editing the user data, updating the role assignment, configuring the messenger access and deleting the delegated user.


The option to delegate control of the account to another user is available by filling in the 3rd party email address of the user to whom the access should be granted and clicking on the **add manager** button:

The option to browse the managed users is available under User profile management menu > **Managed users**:




Managed Users						
<input type="text"/>						<input type="text" value="Search"/>
NAME ↑	E-MAIL ↓	ORGANIZATIONS	MANAGED BY	TOKENCREATED ↑	MESSENGER	
Dr. J. K. Doe	dr.j.k.doe@gmail.com	MDR_test, Testing	john.doe@gmail.com	✓ Sep 6, 2023, 11	<input type="checkbox"/>	...

The option to update the managed user presence in the organizations is available by clicking on the available organizations and switching the correspondent organization switch:




The option to cancel the delegation of the user is available by clicking on the managed by user > Disconnect the manager from the user > Confirm button:



12.1.4 User interface language

The option to change the User interface language is available under User profile management > Languages > select language:




The following languages are available:

- English;
- Ukrainian;
- Russian.

12.1.5 Hiding sensitive information

The option to hide sensitive information (the **patient's** and **uploader names**, **ECG file name** under Tasks section) is available under User profile management > **I'm in public place** switch:



When enabled, all the sensitive information will be blurred during the active session.


12.2 Organizations overview

The Organization section enables a user to enter the organization in order to perform work with patient data inputting and processing.

The following information about the organizations is available under the Organizations section:

Setting	Description
General section	
Name	Indicates the name of the Organization
Path	Indicates the path to the Organization available under the URL of the organization.
Statistics	
Not approved	Indicates the number of Not approved (Pre approved) tasks under the organization.
Processed	Indicates the number of the Approved tasks under the organization.
Cancelled	Indicates the number of the canceled tasks under the organization.
Spending	
Total	Indicates the number of costs the client has spent during working with the software within the organization.



The list of available organizations to user is shown under Organizations screen by the software.



NAME		STATISTICS			SPENDING	
		Not Approved	Processed	Cancelled	Total	
★	MDR_test	1	1	-	-	
	Testing	4	5	-	-	


The option to access the organization is enabled by clicking the Organization:






12.2.1 Organizations filters

Under **Organizations**, the user is enabled to filter the organizations by date:




The user is enabled to set up filter by manually entering the date in DD/MM/YYYY format or via the calendar feature:




12.2.2 Calculate spendings of the organization

The option to calculate spendings of the organization is available by clicking the **Calculate spendings** button:



The option to calculate spendings within all the available organizations is available under **Load Spendings** button:



Spendings calculation depends on the Date **filter**. By default, the spendings are calculated from the first day of the current month.

12.2.3 Export spendings into CSV

The option to export spendings into CSV format is available by clicking the **Export spendings into CSV** button:



Organizations						
<div> <div></div> <div>Load Spendings (31)</div> </div>						
NAME	PATH	STATISTICS			SPENDINGS Total	<div>Export spendings into CSV</div>
		Not Approved	Processed	Cancelled		
★ Testing	testing	4	4	-	-	<div> <div>\$</div> <div> <div></div> <div></div> <div></div> </div> <div></div> <div></div> <div></div> </div>

Note. The option to export spendings becomes available **only** after calculating the spendings. Spendings calculation depends on the Date **filter**. By default, the spendings are calculated from the first day of the current month.

12.2.4 Edit organization

The option to update the organization details is available under the **Edit organization** button:



The following information about the organizations is available to be updated under the Edit organization option:

Setting	Description
General section	
Name	Indicates the name of the Organization
Path	Indicates the path to the Organization available under the URL of the organization.
Visibility level	Indicates the visibility status of the organization to the users within the software. The following visibility levels available: <ul style="list-style-type: none"> Public: sets up the path to “pub_*pathname*” and makes the organization available to operate with, without the user authorization. Private: makes the organization available to operate with only within the users assigned to the organization.
Address	Indicates the physical address of the organization.
Description	Indicates the description of the organization.
Report configuration	
Presets list	Indicates the preset configuration of the report, generated during the ECG task reviewing. The default value of the preset is default . The user is enabled to add, edit and delete the presets.
UTC offset (minutes)	Indicates the timezone of the organization. The option to set the timezone is available by selecting the timezone from the list: <div> EET +03:00 Eastern European Time - Chisinau, Tiraspol, Bați, Bender EET +03:00 Eastern European Time - East Jerusalem, Gaza, Khān Yūni... EET +03:00 Eastern European Time - Helsinki, Espoo, Tampere, Oulu EET +03:00 Eastern European Time - Kyiv, Kharkiv, Odesa, Dnipro EET +03:00 Eastern European Time - Mariehamn EET +03:00 Eastern European Time - Nicosia, Limassol, Larnaca, Stróv... EET +03:00 Eastern European Time - Riga, Daugavpils, Liepāja, Jelgava... </div> By default, the timezone of the organization is EET +03:00 Eastern European Time
Logo	



Logo	Indicates the logo of the organization. Enables a user to load logo if available, change and remove the existing logo of the organization if needed. The following image formats are applicable: .svg, .png, jpeg, .jpg.
------	--

12.2.5 Organization report preset configuration

The option to access organization report preset settings is available under the **Edit organization > Report configuration** section:

Load Spendings (31)

SPENDINGS	Edit organization		
Total			
-	\$	⚙	×

Update organization

Name

Testing

Path *

testing

Visibility level

Private

Address

Test

Description

E-mail: test@xoreserach.com

Report configuration

Presets List

default

Add


Edit

Delete

UTC offset (minutes)

0

Load logo (optional)



×

Cancel

Save

The enabled preset for the organization is set under the **Presets list** dropdown:

Report configuration

Presets List

default

Add

Edit

Delete

UTC offset (minutes)

0

The default value is **default**.

The option to add preset is available under the **Add** button. XOresearch Cardio.AI™ shows the following screen when successful:

Create Report Preset

Name *

Title

Language

English

Time format

HH:mm:ss

20:03:05

Date format

dd MMM

09 Nov

Severity

PQ data

QRS data

QT(c) data

Ordered sections:

Condensed summary

Summary table

Narrative summary

Comments

Daily BPM

Days

Heart Rate Variability (sinus)

ST-segment

Patient's Diary Index

Patient's Diary

Strip Index


Strips

Cancel

Confirm

The following settings are present under **Add** preset and **Edit** preset sections:



Setting	Description
General section	
Name	Indicates the name of the preset.
Title	<p>Indicates the title of the report preset. When the Title is set, it appears on the first page of the task report:</p> 
Language	<p>Indicates the language of the preset being composed on. The following languages are available:</p> <ul style="list-style-type: none"> • English; • Ukrainian; • Russian.
Time format	<p>Indicates the time format of the preset. The following formats are available:</p> <ul style="list-style-type: none"> • HH:MM:SS; • H:MM:SS a.m. / p.m.
Date format	<p>Indicates the date format of the preset. The following formats are available:</p> <ul style="list-style-type: none"> • DD MMM (e.g., 06 Nov); • MMM-DD (e.g., Nov-06); • DD MMMM (e.g., 06 November).
Priority	<p>Indicates the priority of the status. The following priorities available:</p> <ul style="list-style-type: none"> • Highest; • High; • Medium; • Low; • Lowest.
PQ data	Enables to show the PQ interval data with the time between the start of the P wave and the start of the QRS complex
QRS data	Enables to show the QRS complex data with the duration and morphology of each QRS complex.
QT(c) data	Enables to show the QT(c) (Corrected QT interval) data with the



	duration of the QT interval adjusted for heart rate variability.
Ordered sections	Indicates the parts of the report. The order of ordered sections is available to be changed.
Condensed summary	Enables the section which indicates the overall findings and key measurements of the ECG monitoring, including heart rate data, the presence of atrial or ventricular tachycardia, and the burden of ectopic beats.
Summary table	Enables the section which provides a comprehensive overview of key ECG metrics, such as heart rate variability, PQ intervals, and QRS complex durations, summarized in tabular form.
Narrative summary	Enables the section which presents a detailed narrative account of the monitoring period, highlighting significant events, rhythm analysis, and any episodes of bradycardia or tachycardia.
Comments	Enables the section which offers specific observations and insights from the analyzing physician regarding ectopic events, conduction blocks, and other notable findings from the ECG data. This section is a free field to enter the comments during reviewing the ECG analysis.
Daily BPM	Enables the section which charts the daily variations in beats per minute, including maximum, average, and minimum heart rates, as well as occurrences of atrial fibrillation or ventricular blocks.
Days	Enables the section which breaks down the ECG data on a day-by-day basis, allowing for detailed examination of heart rate patterns and ectopic beat occurrences across different times.
Heart Rate Variability (sinus)	Enables the section which displays measures of heart rate variability, offering insights into the autonomic regulation of heart rate during the monitoring period.
ST-segment	Enables the section which visualizes ST-segment deviations and provides analysis on potential ischemic events or abnormalities detected throughout the monitoring duration.
Patient's Diary Index	Enables the section which indexes significant events or symptoms reported by the patient in the diary, correlating them with ECG findings for contextual analysis
Patient's Diary	Enables the section which contains entries from the patient regarding symptoms, activities, or any noteworthy events that may correlate with the ECG data analysis.
Strip Index	Enables the section which organizes the ECG strip recordings by time

	and type of event, facilitating quick access to specific segments of interest for detailed review.
Strips	Enables the section which presents the actual ECG strips that highlight significant cardiac events or intervals of interest identified during the monitoring period.

The option to add a preset is available by filling the **Name** field and clicking the **Confirm** button.

The option to edit preset is available by selecting the preset under the dropdown, clicking the **Edit** button, adding necessary changes and clicking the **Confirm** button.


The option to remove the preset is available by selecting the preset under the dropdown, clicking the **Delete** button and clicking the **Confirm** button.

Note. **Default** preset cannot be deleted.

12.2.6 Remove organization

The option to remove the organization is available under the **Remove organization > Confirm** button:





12.3 Awaiting For Record section overview

Awaiting For Record section enables a user to create tasks that await ECG records before receiving an ECG signal from the device and being transferred to the Tasks section.

Awaiting For Record in Testing

Bind device to patient

Created by: Created

Enter a date range

Filter

ACTION	STATE	DEVICE ID	RECORDING START	FIRST AND LAST NA	BIRTHDAY	INDICATIONS	TAGS	MESSENGER	CREATED	UPDATED	DURATION
Dec 19, 2024											
File	Unknown	test	Dec 19, 2024, 16:10	Test3	Jan 1, 2000	Test3	*	*	Dec 19, 2024, 17:10	Unbound	
File	Uploading	auonm1	Dec 19, 2024, 09:28	Unknown	-		*	*	Dec 19, 2024, 08:28	Unbound	
Oct 25, 2024											
File	Unknown	tryui	Oct 25, 2024, 14:58	Oct 3, 2024			*	*	Oct 25, 2024, 13:58	Unbound	
File	Unknown	tryui	Oct 25, 2024, 14:58	Oct 3, 2024			*	*	Oct 25, 2024, 13:58	Unbound	
File	Unknown	fdddd	Oct 25, 2024, 14:57	Oct 4, 2024			*	*	Oct 25, 2024, 13:58	Unbound	
File	Unknown	tryui	Oct 25, 2024, 14:56	Oct 10, 2024			*	*	Oct 25, 2024, 13:57	Unbound	
File	Unknown	Truei	Oct 25, 2024, 13:10	Oct 10, 2024		вылечено	*	*	Oct 25, 2024, 12:10	7d	
File	Unknown	ytyyy	Oct 25, 2024, 13:09	Oct 3, 2024		mynd,asdfknik	*	*	Oct 25, 2024, 12:09	7d	
Oct 23, 2024											
File	Unknown	ldprt	Oct 23, 2024, 16:15	May 1, 1981			*	*	Oct 23, 2024, 15:16	3d	

Under **Awaiting For Record** section, the following information is available for the user:

- the available action to operate with tasks.

The following actions available:

File action enables a user to select the ECG data file to upload manually into the system.


Start action enables a user to initiate the ECG collection process by the ECG recorder.



Stop action enables a user to stop the ECG collection process by the ECG recorder and start uploading the data to the system;

- State of the task - indicates the status of the task. The following statuses available:
 - Unknown** - indicates that the system does not see a connected and activated device;
 - Not connected** - indicates that the system has identified an activated device, but it is not connected to the system;
 - Waiting for start** - indicates that the system has identified an activated and connected device and is waiting for a start command (automatic start in 10 seconds);
 - Recording** - indicates that the system has identified a connected device that is recording an ECG;
 - Recording is complete** - indicates that the system has identified a connected device that is recording an ECG;
 - Uploading** - indicates that the system has identified a connected device that is transmitting data to the system;
 - Uploading error** - indicates that an error occurred while downloading data from the device. Check the connection of the reading device;
- Device ID - indicates the ID of the ECG device the uploading is being sent from;
- Recording Start - indicates the date and time when the ECG recording was started;
- First and Last Name - indicates the first and last name of the patient;
- Birthday - indicates the date of birth of the patient;
- Indications - indicates the medical indications for the patient;
- Tags - indicates the tags of the task (e.g., test), available to be found by task tag filter;
- Messenger - enables a user to setup a connection between the task and the user's account in Telegram Messenger;
- Created by - indicates the name of the user created the task;
- Updated - indicates the date and time of task last update;
- Duration - indicates the duration of the task. **Unbound** - indicates that the task has no duration.

The user is enabled to filter the tasks under **Awaiting For Record**. The filters are available above the tasks:



The screenshot shows a header bar for 'Awaiting For Record in Testing'. Below the header, there is a row of filter controls. On the left is a green button labeled 'Bind device to patient'. To its right is a 'Created by:' dropdown menu currently set to 'Created'. Next is an 'Enter a date range' input field with a calendar icon and a close 'X' button. Finally, there is a 'Filter' input field with a close 'X' button. The entire filter row is enclosed in a green border.

The following filters are available:



- Created by: available under **Created** dropdown with the available users via the organization.
- Date range: available under **Enter date range** calendar input field, with the option to enter dates manually or via calendar menu:
- First/Second name - available under **Filter** input field, with the option to enter First/Second name of the user in the field.
- Event - available under **Filter** input field, with the option to enter the Event in the field, starting with @ symbol.



- Tag - available under **Filter** input field, with the option to enter the Event in the field, starting with # symbol;

12.3.1 Messenger access

The option to setup Messenger access is available by clicking the **Plus** button under Messenger column:



The system shows the following screen when successful:

Create Messenger Access

Messenger

Telegram

Link

Close
Generate

The Messenger available for creation access is **Telegram**.

The option to create access is available by selecting **Telegram** under messenger dropdown and clicking the **Generate** button:

Create Messenger Access

Messenger *

Telegram ▼

Link

[Redacted Link]

[Copy Icon] [Share Icon]


[QR Code]

Close Generate


The user is enabled to copy link under **copy link** button.



The user is enabled to access link under **Go to link** button:



The user is obliged to click **Start** button under Telegram Messenger.
Telegram Messenger shows the following screen when accessing link:




The option to add a new event to record via Telegram Messenger is available under **Add event** button. The following events available to add:

- Nothing / Accidental tap: /event_none
- Anxiousness: /event_anxious
- Discomfort / Chest pain: /event_chest_discomfort
- Dizziness: /event_dizziness
- Heart palpitations: /event_heart_racing
- Heart Fluttering: /event_heart_fluttering
- Feeling of increased heartbeat: /event_palpitation
- Premonition of fainting / Weakness: /event_feeling_faint
- Shortness of breath / Shortness of breath: /event_short_of_breath

- Fatigue / Exhaustion: /event_tired
- Other: /event_other

12.3.2 Binding device to patient

The option to bind the ECG recorder device to patient and create a task under **Awaiting For Record** section is available under **Bind device to patient** button:



Bind device to patient

Update user data	Advanced settings
<div>First name</div> <div>Last name</div>	<div>Patient Id</div>
<div> <div>Birthdate</div> <div>Age</div> </div> <div> <div>Gender</div> </div>	<div>Ordering Physician</div>
<div>Indications</div>	<div>Ordering Organization</div>
<div>Ordering Phone</div>	<div>Ordering Address</div>
<div>Device Id *</div>	<div>Report Region</div>
<div>Duration</div>	<div>Device Manufacturer</div>
<div>Unbound</div>	<div>Device Name</div>
<div>Recording start *</div> <div>7 Jan 2025, 17:56:13</div>	<div>Service Name</div>
<div>Recording end</div>	
<div>Time zone</div> <div>EET +02:00 Eastern European Time - Kyiv, Kharkiv, Odesa, Dnipro</div>	
<div>Presets List</div> <div>default</div> <div>Show/Edit</div>	
<div>Notification Email</div>	
<div>Assigned to</div>	
<div>Status</div>	
<div>Cancel</div>	<div>Save</div>

Hide Advanced Settings

Setting	Description
General section	
First name	Indicates the first name of the patient.
Last name	Indicates the last name of the patient.
Birthday	Indicates the date of birthday of the patient in the DD MMM YYYY format. The user is enabled to select the date of birthday under the Calendar view
Age	Indicates the age of the patient. This field is being modified by the system according to the changes with Birthday data.
Gender	Indicates the gender of the patient. The following genders available: <ul style="list-style-type: none"> • female; • male; • undifferentiated.
Indications	Provides the indications of the patient.
Device Id	Indicates the Device ID of the patient from which the ECG data was obtained.
Recording start	Indicates the date and time of the ECG recording start.
Duration	Indicates the duration of the ECG recording. The following values are available: <ul style="list-style-type: none"> • Unbound; • 1d; • 2d; • 3d; • 5d; • 7d; • 14d.
Presets list	Indicates the preset configuration of the report of the task, generated during the ECG task reviewing. The default value of the preset is default . The user is enabled to edit the presets.
Notification Email	Indicates the email address the generated task report will be sent to. The value by default corresponds to the value under User data .
Advanced settings	
Assigned to	Indicates the user the patient is assigned to. The available users



	correspond to the users within the organization.
Status	Indicates the status of the task. The following statuses available: <ul style="list-style-type: none"> • Open; • In progress; • Canceled; • Done.
UTC offset (minutes)	Indicates the timezone of the task. The option to set the timezone is available by filling in the field with the time of the time zone other than Greenwich in minutes. For a time zone west of Greenwich, a minus should be set in front of the number. Example: CET - 120.
Patient Id	Indicates the ID of the patient.
Ordering Organization	Indicates the name of the ordering organization of the patient.
Ordering Physician	Indicates the name of the ordering physician of the patient.

Ordering Phone	Indicates the name of the ordering phone of the patient.
Ordering Address	Indicates the address of the ordering organization of the patient.
Report Region	Indicates the region of the report of the task. The following regions available: <ul style="list-style-type: none"> • US; • Canada; • EU; • Ukraine; • Unknown region.
Device Manufacturer	Indicates the manufacturer of the device from which the ECG data was obtained. The following manufacturers available: <ul style="list-style-type: none"> • Life Signals; • Myant; • Cortrium; • Unknown Manufacturer.
Device Name	Indicates the name of the device from which the ECG data was obtained.
Service Name	Indicates the name of the service of the patient.

After filling the **Device ID** parameter and clicking **Save** button, the task appears under **Awaiting For Record**:





When the ECG recorder device connects to the XOresearch Cardio.AI™, the status of the task switches to **Waiting for start**.

If the ECG recorder was terminated early, press the **Stop** button to end the recording. The recording status will change to **Recording completed**.

The ECG recorder would start the data transmission process to XOresearch Cardio.AI™. The status of the record in the system should change to **Uploading**.

12.3.3 Patient data editing


The option to edit Patient's data within the task is available under **Edit personal data** button.

12.4 Tasks section overview

Tasks section enables a user to observe, edit, add and delete tasks, users and roles within the organization.

The Tasks sections is consists of the following subsections:

- Reviewing - enables the user to operate available tasks;
- Uploading - enables the user to access the advanced options of ECG uploading.



12.4.1 Reviewing subsection overview

Under **Reviewing**, the following information is available for the user:

- the available action to operate with tasks. The following actions are:
 - - indicates the inability to operate with a task due to technical difficulties.
 - **Review** - enables a user to edit the ECG task.
 - **View** - enables a user to observe the ECG task.
 - **PDF** - enables a user to download the report of the ECG task.



ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
Aug 23, 2024							
Review	0m	✓	Open	Arfus	+	Aug 23, 2024, 23:07	1d 00h 00m 00s

- the expiration time of the task - indicates how many business hours is left for the task to expire. By default, 7 business hours are set for user to process the task.

ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
Aug 23, 2024							
Review	0m	✓	Open	Arfus	+	Aug 23, 2024, 23:07	1d 00h 00m 00s

- task priority. Task priority serves a hint for a physician considering the prioritizing of ECG processing. In case the software intelligence detects important abnormalities, it sets the higher priority. The following priorities available: Highest, High, Medium, Low, Lowest, Unknown

ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
Aug 23, 2024							
Review	0m	✓	Open	A	+	Aug 23, 2024, 23:07	1d 00h 00m 00s

- task status. The following statuses available:

Open - indicates that the task is available to be edited and no editing actions were applied.

In progress - indicates that the task is currently in the editing process. The status appears after saving changes to the task editing.

Pre-approved - indicates that the ECG task is pre-approved and is available for further editing.

Done - indicates that the ECG report of the task is available to be downloaded and appears after approving the task.

Canceled - indicates that the ECG task is canceled and is unavailable to be processed.

Error - indicates that the error appeared during ECG task processing after uploading.

ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
Aug 23, 2024							
Review	0m	✓	Open	A	+	Aug 23, 2024, 23:07	1d 00h 00m 00s



- patient's name,

ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
						Aug 23, 2024	
Review	0m		Open	A		Aug 23, 2024, 23:07	1d 00h 00m 00s

- task tags - indicates the tags of the task (e.g., test), available to be found by task tag filter,

ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
						Aug 23, 2024	
Review	0m		Open	A	another tes...	Aug 23, 2024, 23:07	1d 00h 00m 00s

- the date of last task update

ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
						Aug 23, 2024	
Review	0m		Open	A	another tes...	Aug 23, 2024, 23:07	1d 00h 00m 00s

- the duration of the record within the task in time format.

ACTION	EXPIRES/QUEUE	PRIORITY	STATUS	PATIENT NAME	TAGS	UPDATED	DURATION
						Aug 23, 2024	
Review	0m		Open	A	another tes...	Aug 23, 2024, 23:07	1d 00h 00m 00s

The user is enabled to filter the tasks under **Reviewing**. The filter are available above the tasks:

Tasks in Testing							REVIEWING	UPLOADING
Upload File	Priority Filters: <input type="text"/>	Priority <input type="text"/>	Status Filters: <input type="text"/>	Status <input type="text"/>	User Filters: <input type="text"/>	User <input type="text"/>	Enter a date range <input type="text"/>	Filter <input type="text"/>

The following filters are available:

- Priority filters: available under **Priority** dropdown with the following priority filters available: Highest, High, Medium, Low, Lowest, Unknown.
- Status Filters: available under **Status** dropdown, with the following status filters available: Open, In progress, Pre-approved, Done, Canceled, Error.
- Assigned to: available under **Assigned** dropdown with the available users via the organization.
- Date range: available under **Enter date range** calendar input field, with the option to enter dates manually or via calendar menu:
- First/Second name - available under **Filter** input field, with the option to enter First/Second name of the user in the field.



- Event - available under **Filter** input field, with the option to enter the Event in the field, starting with @ symbol.
- Tag - available under **Filter** input field, with the option to enter the Event in the field, starting with # symbol;
- Channel - available under **Filter** input field, with the option to enter the Event in the field, starting with \$ symbol;

12.4.2 Reviewing subsection editing

12.4.2.1 Patient data editing

The user is enabled to edit personal data of the patient created with the task under the **Edit personal data** button:

Editing: KHLQJTJGTFIAG8Y7.edf

Id: 4294968455, Id (HEX): 000000100000487

Update user data

First name Last name

Birthday Age Gender

Indications

Device Id Recording start 1 Jan 1970, 00:00:00

Duration Unbound

Presets List default Show/Edit

Assigned to Status Open

Advanced settings

UTC offset (minutes) Patient Id

Ordering Physician Denis Test

Ordering Organization

Ordering Phone Ordering Address

Report Region Device Manufacturer

Device Name Service Name

Hide Advanced Settings

Cancel Save

Editing menu consists of the **Update user data** section and **Advanced settings** section.

The following settings available to be edited under the **Editing** menu:

Setting	Description
General section	



First name	Indicates the first name of the patient.
Last name	Indicates the last name of the patient.
Birthday	Indicates the date of birthday of the patient in the DD MMM YYYY format. The user is enabled to select the date of birthday under the Calendar view
Age	Indicates the age of the patient. This field is being modified by the system according to the changes with Birthday data.
Gender	Indicates the gender of the patient. The following genders available: <ul style="list-style-type: none"> • female; • male; • undifferentiated.
Indications	Provides the indications of the patient.
Device Id	Indicates the Device ID of the patient from which the ECG data was obtained.
Recording start	Indicates the date and time of the ECG recording start.
Duration	Indicates the duration of the ECG recording. The following values are available: <ul style="list-style-type: none"> • Unbound; • 1d; • 2d; • 3d; • 5d; • 7d; • 14d.
Presets list	Indicates the preset configuration of the report of the task, generated during the ECG task reviewing. The default value of the preset is default . The user is enabled to edit the presets.
Notification Email	Indicates the email address the generated task report will be sent to. The value by default corresponds to the value under User data.
Advanced settings	
Assigned to	Indicates the user the patient is assigned to. The available users correspond to the users within the organization.
Status	Indicates the status of the task. The following statuses available: <ul style="list-style-type: none"> • Open;



	<ul style="list-style-type: none"> • In progress; • Canceled; • Done.
UTC offset (minutes)	Indicates the timezone of the task. The option to set the timezone is available by filling in the field with the time of the time zone other than Greenwich in minutes. For a time zone west of Greenwich, a minus should be set in front of the number. Example: CET - 120.
Patient Id	Indicates the ID of the patient.
Ordering Organization	Indicates the name of the ordering organization of the patient.
Ordering Physician	Indicates the name of the ordering physician of the patient.
Ordering Phone	Indicates the name of the ordering phone of the patient.
Ordering Address	Indicates the address of the ordering organization of the patient.
Report Region	Indicates the region of the report of the task. The following regions available: <ul style="list-style-type: none"> • US; • Canada; • EU; • Ukraine; • Unknown region.
Device Manufacturer	Indicates the manufacturer of the device from which the ECG data was obtained. The following manufacturers available: <ul style="list-style-type: none"> • Life Signals; • Myant; • Cortrium; • Unknown Manufacturer.
Device Name	Indicates the name of the device from which the ECG data was obtained.
Service Name	Indicates the name of the service of the patient.

12.4.2.2 Channels editing

The option to edit channels is available under the **Edit channels** button:





Choose leads configuration preset

Some preset name ▼ Delete

Insert preset name Save

Channel 1

Choose lead name

MDC_ECG_LEAD_ES ▼ Invert

Channel 2

Choose lead name

MDC_ECG_LEAD_AS ▼ Invert

Channel 3

Choose lead name

MDC_ECG_LEAD_AI ▼ Invert

Cancel Save

The visibility of the available channels depends on the method of ECG recording and the signal setting.

The following information can be changed under the **Editing channels** menu:



- Preset name of the leads configuration:

Some preset name ▼

Delete

- Proposed preset name field;

Insert preset name

Save

- Name of the channel(s):

Channel 1

Choose lead name

MDC_ECG_LEAD_V6

Invert



Channel 2

Choose lead name

MDC ECG LEAD A

Invert




Channel 3

Choose lead name

MDC ECG LEAD D

Invert



- Inverting the signal of the channel:

Channel 1

Choose lead name
MDC_ECG_LEAD_V6

Invert

Channel 2

Choose lead name
MDC ECG LEAD A

Invert

Channel 3

Choose lead name
MDC ECG LEAD D

Invert

The following lead (channel) names available:

- MDC_ECG_LEAD_I;
- MDC_ECG_LEAD_II;
- MDC_ECG_LEAD_III;
- MDC_ECG_LEAD_AVR;
- MDC_ECG_LEAD_AVL;
- MDC_ECG_LEAD_AVF;
- MDC_ECG_LEAD_V1;
- MDC_ECG_LEAD_V2;
- MDC_ECG_LEAD_V3;
- MDC_ECG_LEAD_V4;
- MDC_ECG_LEAD_V5;
- MDC_ECG_LEAD_V6;
- MDC_ECG_LEAD_ES;
- MDC_ECG_LEAD_AS;
- MDC_ECG_LEAD_AI;
- MDC_ECG_LEAD_A;
- MDC_ECG_LEAD_D.



The option to save preset is available by filling the **Preset name** field, making changes and clicking the upper **Save** button:

Choose leads configuration preset

Some preset name ▼

Delete

Insert preset name

Test

Save

The option to apply changes to the task is available after clicking the lower **Save** button:

Channel 3

Choose lead name

MDC ECG LEAD A ▼


Invert

ECG waveform display

Cancel Save

12.4.2.3 Task reclassification


The option to reclassify task is available under task **options** > **Redo classification (overwrites data!)** button > **Confirm** button:



Note. The reclassification process will overwrite the existing data of the task (e.g., settled annotations)

12.4.2.3 Replacing task data


The option to reclassify task is available under task **options** > **Replace task data (overwrites data!)** button > select ECG file:



Note. The replacing process will overwrite the existing data of the task (e.g., settled annotations)

12.4.2.4 Task cancellation


The option to cancel task is available under task **options** > **Cancel task** button > **Confirm** button:



Note. The task cannot be edited after cancellation. The option to revert cancellation is available by **Reclassification** the task.

12.4.2.5 Task deletion

The option to delete task is available under task **options** > **Delete task** button > **Confirm** button:



12.4.3 Uploading subsection overview

Uploading subsection is showing the uploadings of the ECG data only if the **Auto-process** function is off:

Cardio|AI

Tasks in MDR_test

Upload File Some preset name: Confirm All Auto process Search

REVIEWING UPLOADING 1

Under **Uploading** the following information is available to user:

- The name of ECG file:

	FILE NAME	FIRST NAME	LAST NAME	ASSIGNED TO	AGE	WEIGHT	HEIGHT
Confirm	KHLQTJGTFIAG8Y7.edf	First name	Last name	Assigned to	Age 0	Weight 0	Height 0

- The First Name of the patient:

	FILE NAME	FIRST NAME	LAST NAME	ASSIGNED TO	AGE	WEIGHT	HEIGHT
Confirm	KHLQTJGTFIAG8Y7.edf	First name	Last name	Assigned to	Age 0	Weight 0	Height 0

This setting may be modified before the confirmation step.

- The Last Name of ECG patient:

FILE NAME	FIRST NAME	LAST NAME	ASSIGNED TO	AGE	WEIGHT	HEIGHT
Confirm KHLQJTJGTFIAG8Y7.edf	First name	Last name	Assigned to	Age 0	Weight 0	Height 0

This setting may be modified before the confirmation step.

- Assigned To data:

FILE NAME	FIRST NAME	LAST NAME	ASSIGNED TO	AGE	WEIGHT	HEIGHT
Confirm KHLQJTJGTFIAG8Y7.edf	First name	Last name	Assigned to	Age 0	Weight 0	Height 0

This setting may be modified before the confirmation step.

- Age data:

FILE NAME	FIRST NAME	LAST NAME	ASSIGNED TO	AGE	WEIGHT	HEIGHT
Confirm 09-13-14.EDF	First name	Last name	Assigned to	Age 0	Weight 0	Height 0

This setting may be modified before the confirmation step.

- Weight:


FILE NAME	FIRST NAME	LAST NAME	ASSIGNED TO	AGE	WEIGHT	HEIGHT
Confirm 09-13-14.EDF	First name	Last name	Assigned to	Age 0	Weight 0	Height 0

This setting may be modified before the confirmation step.


The user is enabled to **Edit channels** of the task under corresponding button:

Tasks in Testing							REVIEWING	UPLOADING 1
Upload File	Some preset name:	Confirm All	<input checked="" type="checkbox"/> Auto process				Search	<input type="text"/>
FILE NAME	FIRST NAME	LAST NAME	ASSIGNED TO	AGE	WEIGHT	HEIGHT	Edit channels	
Confirm 09-13-14.EDF	First name	Last name	Assigned to	Age 0	Weight 0	Height 0	Settings Delete	

The user is enabled to **Delete task** under corresponding button > **Confirm**:



User is enabled to search for the tasks under **Uploading** section by using the **Search** field with **File Name** criteria:




12.5 Users section

12.5.1 Users section overview

Users section enables a user to create, invite, manage and delete a user within the organization.


The option to access Users is available under the **Users** tab within the organization:



The following settings are available under the **Users**:

Setting	Description
Active	Indicates the user's activation status switch. When active, the user is functioning in the organization.
User Name	Indicates the name of the user.
E-mail	Indicates the email of the user.
Role	Indicates the role of the user. The available roles of the organization correspond to the roles under the Roles section. The default available roles are: <ul style="list-style-type: none"> • Uploader; • Editor; • Admin.
Created	Indicates the date and time of creation of the user.

The user is enabled to filter the users' data under the **Users** section by the roles, under the **Filter groups** drop-down:



The available filter options correspond to the roles users assigned to.



12.5.2 User creation

The option to create a user within the organization is available under the **Create user** button:

Users in MDR_test

Create user

Invite user

Filter groups

ACTIVE

USER NAME

E-MAIL

ROLE

CREATED

XOresearch Cardio.AI™ shows the following screen when successful accessing:

Create user

First and Last name* *

Required field

Email *

Password *

Select role *

Company name

Contact phone

Contact address

Managed by

Active

Cancel

Create



Setting	Description
First and Last Name	Enables to set the First and Last Name of the user. This field is required .
Email	Enables to set the email of the user. This field is required .
Password	Enables to set the password of the user. The password must include a minimum of 8 characters, comprising special characters, numbers, uppercase letters, and lowercase letters. This field is required .
Select role	Enables to set the role of the user. The available roles correspond to roles under the Roles section. The default roles are the following: <ul style="list-style-type: none"> • Uploader; • Editor; • Admin. This field is required .
Company name	Enables to set the name of the company of the user.
Contact phone	Enables to set the number of the contact phone of the user.
Contact address	Enables to set the address of the user.
Managed by	Enables to set the manager of the user. The available managers correspond to the users within the organization.
Active	Enables to activate or deactivate the user.

The option to create users is available by filling the required fields and clicking the **Create** button:

Create user

First and Last name* *

Test

Email *

test@cardio.ai

Password *

.....

Select role *

Uploader

Company name

Contact phone

Contact address

Managed by

Active

Cancel

Create

12.5.3 User invitation


XOresearch Cardio.AI™ enables the user to invite the user previously created in the system to the current organization. The user is enabled to invite the user by clicking the **Invite user** button > enter email of the user and select the role > **Invite** button:





12.5.4 User editing

The option to edit the user is available under the **Edit user** button:



User editing settings correspond to the User creation settings.

Note. The option to setup User management by another user is not available when editing the user.

12.5.5 User role assignment deletion

The option to remove the user from the organization is available by removing the role assignment of the user from the organization. The option to delete the role assignment is available under **Delete role assignment > Confirm** button:

Users in MDR_test

Create user

Invite user

Filter groups

ACTIVE	USER NAME	E-MAIL	ROLE	CREATED	
<input checked="" type="checkbox"/>	Denis Test		ECG Editor	Apr 19, 2024, 16:00	<div>Delete Role Assignment</div> <div></div>

Do you confirm the deletion of the role assignment of
 Denis Test ?

Cancel

Confirm

12.6 Roles section

12.6.1 Roles section overview

Users section enables a user to create, manage and delete a role within the organization. The option to access Roles section is available under the **Roles** tab within the organization:

Organizations

Tasks

Awaiting For Record

Users

Roles

Roles in MDR_test

Create role



Role Name	Members	Dashbo...					Dashbo...			Manag...			
		View O...	View AL...	Upload...	Edit Ta...	Change...	Change...	ECG Vie...	ECG Re...	Report ...	Organi...	Users ...	Roles ...
admin	2	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
ECG Editor	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Uploader	0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

XOresearch Cardio.AI™ creates a predefined set of roles while creating the organization. The default created roles are: Admin, ECG Editor and Uploader.




Roles section dashboard contains the following components:





Setting	Description
Role Name	Indicates the name of the role.
Members	Indicates the number of users with the correspondent role.
Dashboard	
View Own Tasks	Enables a user to view the tasks the user uploaded ECG of,





	<p>under the Reviewing subsection of Tasks.</p> 
View All Tasks	 <p>Enables a user to view the tasks initiated by all users within the organization of under the Reviewing subsection of Tasks.</p>
Edit Tasks List	<p>Enables a user to Edit personal data of the patient, Reclassify and Reupload the task data.</p>
Change Reviewer	<p>Enables a user to change Assigned to user of the task. under the Edit personal data of the patient.</p>




	<p>Update user data</p> <div> <div>First name</div> <div>Last name</div> </div> <div> <div> Birthday  </div> <div> Age 0 </div> <div> Gender ▼ </div> </div> <div> Indications </div> <hr/> <div> <div>Device Id</div> <div> Recording start 1 Jan 1970, 00:00:00   </div> </div> <div> Duration Unbound ▼ </div> <div> Presets List default ▼ </div> <div> Show/Edit </div> <div> Assigned to Denys Zakhliebaiev ▼ </div> <div> Status Open ▼ </div>
Change Task Status	<p>Enables a user to change Status of the task under the Edit personal data of the patient.</p>

	<p>Update user data</p> <div> <div>First name</div> <div>Last name</div> </div> <div> <div> <div>Birthdate</div> <div></div> </div> <div> <div>Age</div> <div>0</div> </div> <div> <div>Gender</div> <div>▼</div> </div> </div> <div> <div>Indications</div> </div> <hr/> <div> <div>Device Id</div> <div> <div>Recording start</div> <div>1 Jan 1970, 00:00:00</div> <div></div> <div>×</div> </div> </div> <div> <div>Duration</div> <div>Unbound</div> <div>▼</div> </div> <div> <div>Presets List</div> <div>default</div> <div>▼</div> <div>Show/Edit</div> </div> <div> <div>Assigned to</div> <div>Denys Zakhliebaiev</div> <div>▼</div> <div> <div>Status</div> <div>Open</div> <div>▼</div> </div> </div>
Dashboard	
ECG View Access	<p>Enables a user to access ECG task, observe the annotations created by the AI and observe the ECG report.</p> <div> <div>ACTION</div> <div>EXPIRES/QUEUE</div> <div>PRIORITY</div> <div>STATUS</div> <div>PATIENT NAME</div> </div> <div> <div>View</div> <div>0m </div> <div>✓</div> <div>Open</div> <div>Unknown</div> </div> <div> <div>View</div> <div>0m </div> <div>✓</div> <div>Open</div> <div>Unknown</div> </div>

	
Report Final Approve	<p>Enables a user to approve the task, making the report downloadable.</p> <hr/> 
Management	
Organization's Management	Enables a user to edit and remove the organization.
Users Management	Enables a user to create, invite, manage and remove the users within the organization.
Roles Management	Enables a user to create, manage and remove the roles within the organization.
Billing Management	Enables a user to calculate spendings within the organization.

12.6.2 Role management

The option to create role is available under **Roles** section > **Create role** button:



Roles in MDR_test						
Role Name	Members	Dashb...	View ...	View ...	Uploa...	Edit T...
admin	2	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
ECG Editor	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Uploader	0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
test role	0	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

XOresearch Cardio.AI™ shows the following screen when successful:

Create role

Role name *

Permissions:

- ☐ View Own Tasks
- ☐ View All Tasks
- ☐ Upload ECG
- ☐ Edit Tasks List
- ☐ Change Reviewer
- ☐ Change Task Status
- ☐ ECG View Access
- ☐ ECG Review Access
- ☐ Report Final Approve
- ☐ Organization's Management
- ☐ Users Management
- ☐ Roles Management
- ☐ Billing Management

Cancel

Save

Role becomes created when setting Role name, switching the necessary permissions and clicking the **Save** button.



The option to edit the role is available under **Roles** > select Role > **Edit role** button:

Roles in MDR_test															
Create role															
Role Name	Members	Dashb...						Dashb...			Mana...				
		View ...	View ...	Uploa...	Edit T...	Chang...	Chang...	ECG Vi...	ECG R...	Repor...	Organ...	Users ...	Roles ...	Billing...	
admin	2	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Edit role
ECG Editor	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/> ×
Uploader	0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
test role	0	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

The option to remove the role is available under **Roles** > select Role > **Remove role** button > Confirm button:

Roles in MDR_test

Create role

Role Name	Members	Dashb...						Dashb...			Mana...				
		View ...	View ...	Uploa...	Edit T...	Chang...	Chang...	ECG Vi...	ECG R...	Repor...	Organ...	Users ...	Roles ...	Billing...	
admin	2	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Remove role
ECG Editor	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	 
Uploader	0	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
test role	0	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Do you confirm the deletion of the role **ECG Editor**?

Cancel

Confirm

12.7 ECG data input

The option to upload pre-recorded ECG is available under **Upload File** button or drag-n-drop. The option to upload several pre-recorded ECG placed in a folder is available under **Upload Folder** button:

< Organizations

Tasks

Awaiting For Record

Tasks in Testing

Upload File

Upload Folder

Priority Filters: Priority

ACTION

EXPIRES/QUEUE

PRIORITY

STATUS

PATIENT NAME








XOresearch Cardio.AI™ shows the following screen when successful:



1. Header section - enables a user to manage ECG viewer options, save and approve the ECG.
2. Side editing bar - enables a user to choose the ECG periods;
3. Detailed ECG data section - enables a user to view and edit ECG;
4. A report section - enables a user to observe, edit and export the ECG report.

12.8.1 ECG Viewer Header

Header section of the ECG viewer contains the following information:

Setting	Description
Logo	Indicates the logo of the organization: 
ECG file name	Indicates the name of the ECG file. 
Duration of the ECG record	Indicates the duration of the ECG record in days, hours and minutes if applicable. 
Status of the task	Indicates the status of the task: 
Update date	Indicates the time and date of last update of task data: 
Task tags	Indicates the tags of the task:  The option to add tag is available by clicking under Add tag button:  Or via clicking the existing tag. Cardio AI shows the following screen when successful:

Options

RESET

Main options

LANGUAGE

English

TIME FORMAT

HH:mm:ss

DATE FORMAT

dd MMM

Previewer options

CHANNEL

II

SHOW

AMPLITUDE

x 1 scale

ROWS NUMBER

5

ROW DURATION, S

60

ROW HEIGHT, PX

60

COLOR CODES

Visualizer options

CHANNEL

1

LEAD

I

SHOW

AMPLITUDE

10mm/mV

CHANNEL

2

LEAD

II

SHOW

AMPLITUDE

10mm/mV

CENTER LINE

RR INTERVAL

SPEED

25mm/s

ANNOTATIONS

COLOR CODES

RR DIFF, %

20


RULER REPEATS



1


Close

Setting	Description
Main options	
Language	Enables to set the language of task viewer. The following languages available: <ul style="list-style-type: none">English;Russian;Ukrainian.
Time format	Enables to set the time format of the task data.
Date format	Enables to set the date format of the task data.
Previewer options	
Channel	Enables to select the channel for configuration. The available channels correspond to ECG recording device.
Show	Enables to show or hide the channel.



Amplitude	<p>Enables to configure scaling of the amplitude. The available scales are:</p> <ul style="list-style-type: none"> • x1 scale; • x2 scale; • x3 scale; • x4 scale. <p>The option to change scaling is available under Plus and Minus buttons.</p>
Rows number	<p>Enables to set the number of rows under Previewer. The number of rows available from 1 to 20. The default value is 5.</p> 
Row duration, s	<p>Enables to set the duration of rows, in seconds. The following values available:</p> <ul style="list-style-type: none"> • 30; • 60; • 90; • 120.
Row height, px	<p>Enables to set the height of rows, in pixels. The following values available:</p> <ul style="list-style-type: none"> • 30; • 40; • 50; • 60; • 70; • 80; • 90; • 100; • 110; • 120.
Color codes	<p>Enables to set the color codes of annotations availability under Previewer.</p>
Visualizer options	
Channel	<p>Enables to select the channel for configuration. The available channels</p>

	correspond to ECG recording device.
Lead	Enables to select the channel lead for configuration. The available leads correspond to ECG recording device.
Show	<p>Enables to show the channels under Visualizer.</p> 
Amplitude	<p>Enables to set the amplitude of the leads under Visualizer. The following amplitudes available:</p> <ul style="list-style-type: none"> • 5 mm/mV; • 10 mm/mV; • 20 mm/mV; • 40 mm/mV; • 80 mm/mV; • 160 mm/mV; <p>The option to change Amplitude is available under Plus and Minus buttons.</p>
Center line	<p>Enables to show the Center line under Visualizer:</p> 
RR interval	Enables to show the time interval between two successive R-waves of the QRS signal under Visualizer:

	
Annotations	Enables to show the text codes of annotations under Visualizer .
Color codes	Enables to show the color codes of annotations under Visualizer .
Speed	Enables to set the speed of the record under Visualizer . The following speed options available: <ul style="list-style-type: none"> • 12.5 mm/s; • 25 mm/s; • 50 mm/s; • 100 mm/s.
RR diff, %	Enables to set the the percentage difference between successive R-R intervals. The following values available from 0 to 100
Ruler repeats	

The option to reset changes is available under **Reset** button.

The option to save changes is available under **Save** button.

12.8.1.3 Edit ECG task

The option to make the Record editable is available under the **Edit** button:



12.8.1.4 Close ECG task

The option to close ECG task and revert to **Tasks** section is available under **Close** button:



Share Options Edit **Close** Pre-approve Approve Save

12.8.1.5 Pre-approve ECG task

The option to pre-approve task is available under the **Pre-approve** button:

Share Options Edit Close **Pre-approve** Approve Save

Note. Task pre-approving is available only after **Saving** the task.

12.8.1.6 Approve ECG task

The option to approve ECG task and download the report in PDF format is available under the **Approve** button:

Close Pre-approve **Approve** Save

Note. Task approving is available only after **Saving** the task.

12.8.1.7 Save ECG task

The option to save changes after editing the ECG task is available under **Save** button:

Share Options Edit Close Pre-approve Approve **Save**

XOresearch Cardio.AI™ shows the following notification when successful:

 Analytics data updated


12.8.2 ECG Viewer Editor

ECG Viewer Bulk Editor consists of the following elements:

1. Side editing bar - contains all and classified beats by annotations;
2. Poincare plot - enables to navigate and select beats;
3. Beats list - enables to bulk select and edit beats




4. Beats cluster panel - enables to compare beats via clusters;
5. Beats cross-annotations list - enables to observe and manage beats with multiple annotations.




12.8.2.1 Side editing bar

Under side editing bar, XResearch Cardio.AI™ collects and shows all the beats, normal beats and found annotations. All beats highlighted with **black**, normal and sinus beats highlighted with **green**, annotations highlighted with other colors.





The user is enabled to click on the highlighted segments, and software will focus on that selected segment, including a detailed view of the selected segment:




Multi selection of episodes is available via SHIFT or CTRL buttons.


12.8.2.2 Poincare plot

Poincare plot enables a user to review, view and navigate all recorded beats, including normal, and annotations (abnormalities)

The option to navigate to beat is enabled by clicking on the beat:



The user has the ability to select multiple bits by drawing an area on the Poincare plot by left-clicking and moving the cursor over the plot:




By default, Poincare plot is shown in the mode RR +1. The option to switch to RR-1 mode is available under the corresponding switch:

All beats


Total 97,900 beats | 4,715 episodes




The option to show only normal beats to the left and to the right is available under **only n-n** switch:



The option to filter beats under Poincaré plot is available by clicking on the following filter elements:




The user is enabled to change the length of filter element by moving the filter elements borders by clicking them and moving the cursor:




Filtered beats are displayed under the Beats list:



The option to reset filter is available under **Reset filter** button:




When filtered and selecting the single beat, the option to revert to a filter view is available under the following button:



12.8.2.3 Beats list


Beats, selected under the **Side editing bar** or **Poincare plot** are displayed under the Beats list:



Setting	Description
*	Indicates the number of the beat within the beats under the Poincare plot.
Event	Indicates the name of the annotation correspondent to the beat.
-,%	Indicates the difference in % between the beat and the beat which is left to the beat.
RR,ms	Indicates the distance in ms between the beat and the beat which is left to the beat.
aBPM	Indicates the average BPM of the beat (calculated for the 6 seconds).
Position	Indicates the position (time) of the beat location on the ECG recording








The user is enabled to filter the settings in ascending and descending order by clicking on the setting in the column:




	Event ↓	%	RR,ms	aBPM	Position
1	> N	+1	886		07:47:18 11 Nov
2	N nn	+2	869		07:47:17
3	N nn		846		07:47:17

The option to select beat is available by clicking on the beat. The following options available to navigate within beats:

-  (Spacebar)- enables to select the next beat;
-  (CTRL + Spacebar) - enables to select the previous beat;
-  - enables to select the next twentieth beat;
-  - enables to select the previous twentieth beat;
-  - enables to select the first beat.

The option to access editing menu of the beats under **Beats list** is available by clicking the right button while selecting the beats, or via **Edit** button:




The left screenshot shows a 'SMART ACTIONS' menu with various beat types: UNK, A, N, E, N, R, A, N, S, WSP, VPCF, JESC, G, UNK, REPLACE, ADD, SOFT, Ventricular, Atrial, Junctional, Block, Pre-Excitation, Sinus, Bundle Branch Block, Noise. The right screenshot shows an 'Actions menu for selected beats' with a green box around the edit icon (a pencil).

12.8.2.4 Smart actions

Smart Actions - the Editing menu, enables a user to manage beats within the **Viewer** and **Visualizer** sections of **ECG** task.

The following sections available under **Smart Actions** menu:

- Remove annotation - indicates the existing annotation, applied to the beat. When clicking - removes the annotation and classifies the beat as normal. Available only for **abnormal** beats.
- History of **Smart actions**, with the recently applied annotations:



Note. Smart actions history differs within the number of selected beats (1-3, 4+).

- **Replace** - when enabled, allows to replace the beat with another annotation;
- **Add** - when enabled, allows to add the annotation to the beat. When using **Add** option, the added annotation replaces the previous annotation, if applicable. ;
- **Soft** - when enabled, allows to **soft add** the annotation to the beat. When **soft adding**, the added annotation does not replace the previous annotation, if applicable.;

The following annotations available under **Smart actions**:

- **Ventricular** - contains the following annotations:
VPC - Ventricular Premature Contraction;




- VPCF - Fusion of Ventricular And Normal Beat;
- VPCI - Ventricular Interpolated Beat;
- RONT - R-On-T Premature Ventricular Beat;
- VESC - Ventricular Escape Beat;
- V2 - Ventricular Couplet;
- V3 - Ventricular Triplet;
- VL - Ventricular Bigeminy;
- VTRG - Ventricular Trigeminy;
- VFIB - Ventricular Fibrillation;
- VFLU - Ventricular Flutter;
- VTDP - Torsades De Pointes Ventricular Tachycardia;
- MOV - Monomorphic Ventricular Tachycardia;
- PLVT - Polymorphic Ventricular Tachycardia;
- VR - Idioventricular (Ventricular Escape) Rhythm;
- VAIR - Accelerated Idioventricular Rhythm;
- **Atrial** - contains the following annotations:
 - APC - Atrial Premature Contraction;
 - ABER - Aberrated Beat;
 - NPW - Non-Conducted P-Wave (Blocked);
 - AESC - Atrial Escape Beat;
 - A2 - Atrial Couplet;
 - A3 - Atrial Triplet;
 - ABI - Atrial Bigeminy;
 - ATRG - Atrial Trigeminy;
 - AFIB - Atrial Fibrillation;
 - AFLU - Atrial Flutter;
 - PAT - Paroxysmal Atrial Tachycardia;
 - MAT - Multifocal Atrial Tachycardia;
 - AAT - Automatic Atrial Tachycardia;
 - AERY - Atrial Ectopic Rhythm;
 - WSP - Wandering Sinus Pacemaker Within The Sinus Node;
 - ARYU - Upper Atrial Rhythm;
 - ARYM - Middle Atrial Rhythm;
 - ARYL - Lower Atrial Rhythm;
- **Junctional** - contains the following annotations:
 - JPC - Junctional (Nodal) Premature Contraction
 - JESC - Junctional (Nodal) Escape Beat
 - J2 - Junctional Couplet;
 - J3 - Junctional Triplet;
 - JB - Junctional Bigeminy;
 - JTRG - Junctional Trigeminy;
 - JT - Junctional Tachycardia;
 - RECP - AV Reciprocating Tachycardia;



RNTR - Reentrant AV Nodal Tachycardia;
 WAP - Wandering Pacemaker From The Sinus Node To (And From) The A-V Node;
 IRYE - AV Junctional (Nodal) Escape Rhythm;
 IRYA - Accelerated Av Junctional (Nodal) Rhythm;

- **Block** - contains the following annotations:
 AV1 - First Degree AV Block;
 AV2I - Second Degree AV Block Type I;
 AV2II - Second Degree AV Block Type II;
 AV3 - Third Degree AV Block;
 AVDI - AV Dissociation With Interference;
 AVDS Isorhythmic AV Dissociation;
 AVDC - Complete AV Dissociation;
 SA2I - Second Degree SA Block Type I;
 SA2II - Second Degree SA Block Type II;
 SA3 - Third Degree SA Block;
 PAUS - Pause;
 AV2 - Second Degree Av Block;
- **Pre-Excitation** - contains the following annotations:
 WPWA - Wolf-Parkinson Type A;
 WPWB - Wolf-Parkinson Type B;
 LGL - Lown-Ganong-Levine Syndrome.
- **Sinus** - contains the ARHY - Sinus Arythmia annotation;
- **Bundle Branch Block** - contains the following annotations:
 BBB - Bundle Branch Block Beat (Unspecified);
 LBB - Left Bundle Branch Block Beat;
 LBBI - Incomplete Left Bundle Branch Block Beat;
 RBB - Right Bundle Branch Block Beat;
 RBBI - Incomplete Right Bundle Branch Block Beat;
 BBLA - Left Anterior Fascicular Block Beat (Common);
 BBLP - Left Posterior Fascicular Block Beat (Rare);
 BBBI - Bifascicular Block Beat;
 BBTI - Trifascicular Block Beat
 BBBL - Bilateral Bundle-Branch Block Beat
 BBU - Intraventricular Conduction Disturbance (Non-Specific Block)
- **Noise** - contains the following annotations:
 UNK - Unclassifiable Beat;
 ZZZ - Noise (No Signal);
 Z - Noise Moderate;
 ZZ - Noise Severe;
 A - Artifact.


Visualizer: The option to remove annotation is available by selecting the beat > **Delete selected annotations for selected beats** button:



12.8.2.5 Beats clusters panel

Under the clusters panel, the user is enabled to select clusters of the channels to be displayed the **Poincare plot**.

The option to select channel is available under the **Lead** dropdown:




The option to select cluster is available by clicking on the available cluster:




The number on cluster indicates the number of beats within the cluster.

The user is enabled to filter clusters by the beat number under the following button:




The user is enabled to multi select the cluster by clicking SHIFT and selecting the clusters.

The option to reset the selection filter of the clusters is enabled under the **Reset filter** button:




12.8.2.6 Beats cross annotations list

Under the **cross annotations list**, beats are grouped by the events (annotations), and the quantity:




A2*	1
A3*	2
APC*	2
JPC+V2	1
MAT*	1
MOVT*	1
N*	2
PLVT*	5
V2*	44
V3*	26
VBI*	1
VPC*	6

The user is enabled to select the events by clicking on them. The option to reset the selection is available under the **Reset filter** button:



1/1	
Reset filter	
A2*	1
A3*	2
APC*	2
JPC+V2	1
MAT*	1
MOVT*	1
N*	2
PLVT*	5
V2*	44
V3*	26
VBI*	1
VPC*	6

XOresearch Cardio.AI™ indicates the cross-annotations with the following indication:




A2*	1
A3*	2
APC*	2
JPC+V2	1
MAT*	1
MOVT*	1
N*	2
PLVT*	5
V2*	44
V3*	26
VBI*	1
VPC*	6

The cross-annotations must be reviewed by the healthcare professional.


12.5.3 ECG Viewer Previewer

Previewer of the ECG viewer displays the area in which multiple hearts beats are included:






When enabled in options, Previewer includes the color coded annotations.
The left side of each row displays the date and time of the recording section:




The user is enabled to navigate via the Previewer by clicking on row area:




12.8.4 ECG Viewer Visualizer

Under ECG Viewer Visualizer, the user is able to observe and manage beats and annotations correspond to Previewer and Editor.




Detailed ECG data section contains the preview section of the ECG divided into parts for each minute of the ECG recording for the whole ECG recording period:






The user is enabled to choose the period by scrolling and choosing the required period. Chosen period is highlighted at the bottom detailed viewer section:




The option to select a beat is available by clicking on the beat



The applicable annotation to the beat set by XOresearch Cardio.AI™ is available at the upper side of the beat, and shows the name by hovering the cursor above:






At the visualizer, there is a line at the top: a parameter at the top of the line indicates the distance between the beat and the one on the left; a parameter at the bottom indicates the average BPM:



MD

CE 0123

Under Visualizer, the user is enabled to measure the AMP by clicking the beat > clicking **AMP** button > Click the left mouse button on the Visualizer beat and swipe the cursor up or down:



The user is applicable to set up several AMP measurements. The option to remove AMP measurements is enabled by clicking the X button under **AMP** button.

Under Visualizer, the user is enabled to measure the speed by clicking the **Speed** button > Click the left mouse button on the Visualizer beat and swipe the cursor left or right:




The user is applicable to set up several speed measurements. The option to remove Speed measurements is enabled by clicking the X button under **Speed** button.

The user is enabled to add a strip of the beat to the report by selecting the **beat > Add strip** button:





The option to navigate to a certain time / sample is available under **Goto** button > Select **Date** and **Time** > Enter the **Sample** number > **Go** button:







Go to Time/Sample

Date/Time
Apr 3, 2024, 10:55:29 AM 

Sample
1345 


Cancel 

The option to view PQRS measures recognised by AI is available by double clicking on the beat:




The option to access Editing menu is available by clicking on the right button at the beat:





Under Editing menu, the user is enabled to select the beats by the rhythm regularity. The following actions available:


- Select by rhythm regularity, left:




- Select by rhythm regularity:




- Select by rhythm regularity, right:




The user is enabled to browse the history of recent actions under the following field:



The option to undo changes is available under **Undo** button:




The option to redo changes is available under **Redo** button:



12.5.5 ECG Viewer bird view

XOresearch Cardio.AI™ enables a user to check and navigate to the events during the recorded ECG day and night periods via bird view:




Time appears above the bird view section when hovering over the selected fragment.

12.9 ECG data report

12.9.1 ECG data report overview

The option to view report of the ECG is available under Overview button:





Report is divided into the following sections:

- Personal data section - contains the following information of patient: date of birth, gender, ordering organization, contacts, ordering physician, signature, device id, lead configuration, interpretation physician, recording time, amount of analyzed time, enrollment time;
- Comments section - contains the optional comments; The option to write comments is available by clicking on the **Comments** field and entering the text:

Comments


Test

- Condensed summary - provides a concise overview of key monitoring data for quick reference. It includes essential information about the monitoring duration, predominant rhythm, heart rate statistics, and major findings, such as ectopic beats, blocks, and tachycardia events. The purpose is to offer a high-level snapshot of the monitoring results for easy comprehension.
- Narrative summary - offers a detailed and chronological account of the monitoring session. It presents a comprehensive analysis of the data, including specific events, their durations, and their timestamps. The purpose is to provide healthcare professionals with a thorough understanding of the patient's cardiac activity during the monitoring period, allowing for more in-depth evaluation and decision-making. It also highlights notable episodes and deviations from normal rhythm and presents relevant metrics and measurements.



- Summary - provides a comprehensive and structured overview of the key findings and metrics derived from the cardiac monitoring session. It serves as a consolidated report that healthcare professionals can reference to quickly assess the patient's cardiac health and identify any notable deviations from the norm.
- Daily BPM - provides a BMP (beats per minutes), including ectopic beats.
- BMP (sinus) - provides a BMP on sinus beats, excluding ectopic beats;
- PQRST (sinus) - provides a information about PQ interval, QRS complex, QT/QTc intervals
- Annotations list - provides details of different annotations, according to the timeline. The explanation of the abbreviations is placed below the list. Every annotation has it's own features.
- Heart Rate Variability (sinus) - provides various aspects of heart rate variability and sinus rhythm. They provide insights into the health of the cardiovascular system and the variability in time between successive heartbeats.
- ST-segment and T-wave type - provides the length and direction of ST segment and determines the type of T-wave.
- Strip Index table - contains information about specific cardiac events, including their labels, notes, associated heart rates, and timestamps;
- Strips section - provide additional details or data related to specific events or conditions mentioned earlier. It includes heart rate measurements (in BPM) and timestamps for each event.
- Patient's diary index table - contains information about specific cardiac events highlighted by the patient, including their labels, notes, associated heart rates, and timestamps;
- Patient's diary strips - provides additional details or data related to specific events highlighted by the patient, or conditions mentioned earlier. It includes heart rate measurements (in BPM) and timestamps for each event.

The option to approve report is available under **Approve** button:



The option to export the report is available under Organization view after approving the report > **Download report** button:



Tasks in Testing

Upload File

Upload Folder

Priority Filters:

Priority

Status Filters:

Status

Assigned to:

Assigned

Enter a date range

Filter

ACTION	PRIORITY	STATUS	PATIENT NAME	FILE	ASSIGNED TO	UPLOAD BY	TAGS	UPDATED	DURATION
<div>View</div> <div>PDF</div>	Done	Unknown (age 55)	A	S	S	Feb 19, 2024, 15:44	23h 59m 50s		
<div>View</div> <div>PDF</div>	Done	V A (age 64)	9	P	S	Feb 19, 2024, 15:43	3d 00h 03m 19s		
<div>Review</div> <div>0m</div>	Open	K E (age 36)	4	Unknown	S	Feb 19, 2024, 15:41	2d 00h 01m 39s		

Note. After task approving, the report is automatically send to the email address stated under **Notification email** setting under **User profile management**:

"Cardio.AI" patient report: Inbox

sender@sapiensapi.com
to me

4:25 PM (0 minutes ago) ☆ 📧 ↻ ⋮

CardioAI[®]
MONITORING

Report is ready

Information about the patient

Full name:

Date of Birth: —

Gender: Unknown

Ordering Organization:

Contacts: ,

Ordering Physician:

Device ID:

Recording Time: 3h 26m

Enrollment Time: 01 Jan 1970 00:00:00 CET - 01 Jan 1970 03:26:39 CET

Performer: Testing

Legal Address: Test

Interpretation Physician:

Risk levels

Priority level	Description
Ok	Your heart activity is normal, it is recommended to repeat the test after a year and a half or sooner.
Low	Minor violations that do not require further consultation or treatment.
Moderate	Identified violations that may require a doctor's attention. It is recommended to consult within a month.
Important	Serious violations. Consultation with a cardiologist is necessary within a week.
Critical	Life-threatening violations. See a doctor immediately!

View the report

Cardio.AI | Innovative solutions for heart health

Visit the site

The user is able to navigate to heartbeat sections of interest by clicking on cells in the report:





label

The user is able to observe the data from the ECG channels under **Strips** section of the report > Labels. The table of labels is available under **Strip Index** section:



Strip Index

Label	Note	BPM	Time
Sinus BPM Max			17 Aug 22:14:30
Sinus BPM Min			18 Aug 10:16:31
Atrial Premature Contraction			17 Aug 21:17:26
Junctional (Nodal) Premature Contraction			18 Aug 09:56:16
Aberrated Beat			18 Aug 10:18:05
Non-Conducted P-Wave (Blocked)			18 Aug 17:45:23
Ventricular Premature Contraction			17 Aug 19:52:02
Junctional (Nodal) Escape Beat			18 Aug 14:11:21
Sinus Arrhythmia			17 Aug 19:59:23
Wandering Sinus Pacemaker Within The Sinus Node			18 Aug 13:25:21
Wandering Sinus Pacemaker Within The Sinus Node			18 Aug 13:26:00
Atrial Ectopic Rhythm			18 Aug 14:51:05
Atrial Ectopic Rhythm			18 Aug 14:52:10
Atrial Bigeminy			18 Aug 17:44:41
Atrial Flutter			18 Aug 11:19:03
Atrial Flutter			18 Aug 11:19:34
Atrial Flutter			18 Aug 11:22:15
Atrial Flutter			18 Aug 11:23:42
AV Junctional (Nodal) Escape Rhythm			18 Aug 09:45:47
First Degree AV Block			18 Aug 10:28:29
Second Degree SA Block Type I			18 Aug 17:31:50
Lown-Ganong-Levine Syndrome			18 Aug 11:15:54
Lown-Ganong-Levine Syndrome			18 Aug 11:16:07
Pause			18 Aug 07:55:47
Atrial Couplet			17 Aug 22:32:36
Atrial Triplet			18 Aug 14:07:40
Nonsustained Atrial Flutter			18 Aug 11:17:49
Nonsustained Atrial Flutter			18 Aug 11:18:18
Nonsustained Atrial Flutter			18 Aug 11:18:29
Nonsustained Atrial Flutter			18 Aug 11:27:14

The user is enabled to navigate to the strip by clicking on strip under **Strip Index** table of labels.

By default, the data under **Labels** is being shown from the **ES, AS, AI** channels.




Strips

Sinus BPM Max

17 Aug 22:14:30

12.5 mm/s 110 BPM




The option to expand the data from all channels is available under **Strip 1** button:







Note. The option to expand the channels is depend on the ECG data source and the availability of channels from the ECG recording device.
By default, the speed is 12.5mm/s. The option to expand the amplitude is available by clicking the entry:



The option to access and share the entire ECG record to observe for the 3rd party person is available under **See the entire ECG** link:


Def uk  Highest		Date of Birth 22 Jan 1997 (21 yrs)	Gender Male	Testing Org Address Description Testing organization111
Ordering Organization Organization		Device ID	Recording Time 23h 59m	Enrollment time 17 Aug 2018 19:43:00 18 Aug 2018 19:42:50
Contacts phone address		Lead Configuration EASI leads	Analyzed (noise skipped) 23h 59m	
Ordering Physician Superuser Signature		Interpretation Physician Superuser Signature	POWERED BY 	

Note. The link is workable within the 90 days since report generation.

12.9.2 ECG data report sections managing

The option to manage ECG data report sections is available under the **Report menu** section:





XResearch Cardio.AI™ shows the following screen when successful:



Update Report Preset

Presets List

default

Title

Language

English

Time format

HH:mm:ss

20:03:05

Date format

dd MMM

09 Nov

Priority

PQ data

QRS data

QT(c) data

Ordered sections:

Condensed summary

Summary table

Narrative summary

Comments

Daily BPM

Days

Heart Rate Variability (sinus)

ST-segment

Patient's Diary Index

Patient's Diary

Strip Index

Strips

Cancel

Confirm

The settings under **Report menu** correspond to the settings under [Report preset configuration](#).

12.9.3 ECG data report editing

XOresearch Cardio.AI™ enables a user to edit the following sections of the report inside the ECG task:

- Condensed summary;
- Narrative summary;
- Comments.

The option to edit the sections above is available by clicking the section, or by clicking the **Edit** button:



Condensed summary

The monitoring was worn from **31 Dec 1969 19:00:00** for **3h 26m**.

The predominant rhythm shows *Normal Sinus Rhythm*.

The findings of the monitor are detailed below:

- The minimum heart rate was 54, the maximum heart rate was 115, and the average heart rate of 73.
- There were 2 PVCs (<0.1% burden).
- There was 2 heart block (<0.1% burden) and 1 significant pauses.

The option to remove the data entry is available by clicking the **Remove** button:

Condensed summary

The monitoring was worn from **31 Dec 1969 19:00:00** for **3h 26m**.

The predominant rhythm shows *Normal Sinus Rhythm*.

The findings of the monitor are detailed below:

- The minimum heart rate was 54, the maximum heart rate was 115, and the average heart rate of 73.
- There were 2 PVCs (<0.1% burden).
- There was 2 heart block (<0.1% burden) and 1 significant pauses.

13. Data Input and Output:

Data Input:

- XOresearch Cardio.AI™ accepts ECG file data in the following formats: EDF, BDF.
- Ensure that all input data is accurate and complete;

Data Output:

- XOresearch Cardio.AI™ generates reports based on analysed ECG data and displays it on the screen on purpose. The user is enabled to export this report as a PDF report for sharing with other healthcare professionals.

14. User Authentication and Access Control:

User Authentication: Each authorized user is required to log in using their unique username and password. It is essential to keep login credentials confidential. Login credentials are being provided by the XOresearch SIA directly, via the contact email, or via the contact webform under XOresearch Cardio.AI™ [website](#).

Access Control: The software offers role-based access control, ensuring that users only have access to the features and patient data relevant to their role. Administrators can manage user permissions.

There are 4 types of users to access XOresearch Cardio.AI: Support, Administrator, Editor, and Uploader. A brief description of each of them is given below.



Support: This is the user responsible for managing organizations (hospitals or clinical settings) and user profiles within these organizations. Only XOresearch personnel can have this type of access.

Uploader: This is a user who can upload ECG data and download the report to be delivered to a patient inside the organization.

ECG Editor: This is a user with Uploader access and a few more permissions.

Administrator: This is the user with an admin role inside a given organization.

User Type	User permissions
Uploader	<ul style="list-style-type: none"> • Upload ECG records; • Create tasks based on uploaded ECG records; • Manage metadata for the created tasks; • View only the created tasks;
ECG Editor	<ul style="list-style-type: none"> • Upload ECG records; • Create and manage tasks based on uploaded ECG records; • View, edit ECG, create, manage and export reports for the ECG tasks within the organization; • Manage metadata for the tasks within the organization.
Admin	<ul style="list-style-type: none"> • Upload ECG records; • Create and manage tasks based on uploaded ECG records; • View, edit ECG, create, manage and export reports for the ECG tasks available within the organization; • Manage metadata for the tasks within the organization; • Manage users, roles and permissions within the organization.
Support	<ul style="list-style-type: none"> • Upload ECG records; • Create and manage tasks based on uploaded ECG records; • View, edit ECG, create, manage and export reports for the ECG tasks available within the organizations; • Manage metadata for the tasks within the organizations;



	<ul style="list-style-type: none"> • Manage users, roles and permissions within the organization; • Managing the organizations, users, roles and permissions within the software.
--	---

Note: allocation of 'customizable' permissions is the responsibility of the health institution's admin.

The Support role is intended to be used only by the XOresearch Cardio.AI™ staff members.

15. Data Security and Privacy:

XOresearch SIA places the utmost importance on the security and privacy of patient data. We employ industry-standard encryption protocols to ensure the confidentiality and integrity of patient data during both transmission and storage. Additionally, our software complies with all relevant data privacy regulations, including but not limited to the Regulation (EU) 2016/679 (General Data Protection Regulation - GDPR) and the Health Insurance Portability and Accountability Act of 1996 (HIPAA). These measures are in place to safeguard patient privacy and data security.

Additional User Security Recommendations:

In addition to the security measures we have implemented, we strongly recommend that users take the following steps to enhance cybersecurity while using XOresearch Cardio.AI™:

Keep Your Login Information Secure: Never share your login credentials, and ensure they remain confidential. Avoid writing down login information or storing it near your computer.

Access Control: Always log out of XOresearch Cardio.AI™ when not actively using it, especially when in shared or public environments.

Regularly Change Your Password: Change your password at the first login and periodically thereafter. Use strong passwords that include a minimum of 8 characters, comprising special characters, numbers, uppercase letters, and lowercase letters.

Avoid Common Passwords: Refrain from using easily guessable passwords, such as simple combinations or common words. Never use the same password for multiple devices or accounts.

Verify Website URLs: Always verify the URL address before logging into any site. Secure websites start with "https," and a green lock symbol should be displayed in the URL bar.

Install Antivirus and Antispyware Software: Protect your computer by installing and regularly updating antivirus and antispyware software.



Report Suspicious Activity: If you notice any unexpected behavior on your system while using XOresearch Cardio.AI™, please contact our support team. If necessary, we will notify you via email and/or our website if the system faces potential threats that require downtime for resolution.

System Updates: Regularly update your browser used to access the XOresearch Cardio.AI™ and any associated systems to apply the latest security patches. This is crucial to safeguard against newly identified vulnerabilities.

Data Consent: Obtain explicit patient consent before storing or processing data with XOresearch Cardio.AI™, especially for long-term storage or data sharing with other entities. Document consent as part of the patient's medical record.

Anonymization Best Practices: For all identifiable patient data, follow anonymization protocols to prevent unauthorized access. This includes restricting access to only authorized personnel and applying anonymization techniques where applicable, especially when data is shared outside the organization.

Continuous Improvement and User Notifications:

As part of our commitment to security, we continuously monitor cybersecurity threats and make necessary improvements. We will keep you informed of software updates, revisions, or additional security measures through email notifications, ensuring that you have access to the latest safeguards and enhancements.

16. Troubleshooting:

If you encounter technical issues or unexpected errors while using XOresearch Cardio.AI™, please contact our technical support team at getintouch@xoresearch.com.

17. Availability of the Instructions for Use (IFU):

The Instructions for Use (IFU) for XOresearch Cardio.AI™ is provided in electronic format.

The electronic version (eIFU) is available for observation from the official SIA XOresearch Support Centre website at: <https://support.cardio.ai/ifu/index.html>. The version number and revision history of the IFU are documented above in this document to ensure full traceability.



Users can request an additional copy by contacting XOresearch Support via email at getintouch@xoresearch.com.

Users can request a paper copy of the Instructions for Use. Requests must be made through the manufacturer's official contact channels (email: getintouch@xoresearch.com | Phone: +371-67-305-084). The IFU will be shipped within 7 calendar days of receiving the request or provided at the time of device delivery if requested during order placement.

It is the responsibility of the user to ensure that they are referring to the latest version of the IFU, which can be verified on the XOresearch website.

The manufacturer maintains an effective update notification system. Users who have accessed the IFU online are advised to check the official website periodically for updates. In case of safety-related updates or corrective actions, registered users will be informed via email notifications.

According to Article 5(8) of Regulation (EU) 2021/2226, XOresearch has implemented a system to clearly indicate when the Instructions for Use have been revised and to inform users if the revision was necessary for safety reasons. This system is maintained in compliance with the QMS document control procedures, which include version tracking, change history, and registered user notification workflows.

18. Limitations

XOresearch Cardio.AI™ is a clinical decision support software designed to assist healthcare professionals in ECG data analysis. While using, the following limitations should be considered:

Clinical Decision Support Only

XOresearch Cardio.AI™ does not provide a definitive diagnosis and is not intended to replace clinical judgment. It serves as an aid to qualified healthcare professionals who must interpret the results in the context of the patient's clinical presentation.

Dependence on Input Data Quality

The accuracy of analysis depends on the quality and integrity of the ECG data. Incorrect lead placement, signal noise, or incomplete recordings may affect performance and lead to misinterpretation.

No Real-Time Monitoring or Emergency Alerts

The software processes ECG data retrospectively and does not support real-time monitoring or automated alerts for critical cardiac events. It is not intended for use in emergency decision-making.



Pacemaker Signal Limitations

The software does not reliably detect or differentiate ECG signals originating from implanted pacemakers or defibrillators. It cannot be used as a tool for patients with these devices.

ECG Format Compatibility

XOresearch Cardio.AI™ supports ECG data import in EDF and BDF formats only. ECG recordings in other proprietary formats may not be compatible unless converted to a supported format.

Regulatory Scope and Intended Use

The software is classified as a Class IIa medical device under MDR (EU) 2017/745 (Rule 11). Its intended use is limited to the scope defined in the regulatory documentation and certification. Any use beyond this scope is not covered by the manufacturer's intended purpose.

System and Environmental Requirements

XOresearch Cardio.AI™ is a web-based application requiring stable internet connectivity and a compatible browser (Google Chrome 116+, Microsoft Edge 126+, or Opera 113+). Performance may be affected if system requirements are not met.

User Training Requirement

The software should only be used by qualified healthcare professionals who have reviewed the Instructions for Use (IFU) and completed appropriate training. Improper use may result in misinterpretation of ECG data.

Risk of False Positives/Negatives

Despite rigorous validation, the software may produce false-positive or false-negative classifications. Clinical verification of AI-generated annotations is **required** before making patient management decisions.

Data Storage and Retention

ECG data is stored for a limited period per the manufacturer's data retention policy. Users must comply with applicable data protection regulations regarding the storage, processing, and transfer of patient information.

System Responsiveness During Large ECG File Uploads

When uploading large ECG data files, the system may appear temporarily unresponsive while preparing the files for upload. This is a known technical limitation due to browser processing constraints. However, this does not affect data integrity or the accuracy of ECG analysis.



To continue working without interruption, users can open a new browser tab while the upload progresses in the background. No data loss or functional disruption occurs, and all uploaded files are processed as expected.

19. Manufacturer's Declaration

We, SIA XOresearch, declare that this Instructions for use accurately represents the use and troubleshooting procedures for XOresearch Cardio.AI™.

Any serious incident related to the device must be reported to XOresearch SIA and to the competent authority of the Member State in which the users and/or patients are established.

